

## SAM target sgRNA cloning protocol – S. Konermann, Zhang lab, 2014

Optimized sgRNAs for any coding human gene can be found using our SAM Cas9 activator design tool:

<http://sam.genome-engineering.org/database/>

In order to clone the guide target sequence into the sgRNA(MS2) cloning backbone (addgene #61424 ) or lenti sgRNA(MS2)\_zeo backbone (addgene #61427), synthesize two oligos of the following form. **Both plasmids have the same overhangs:**



*Example oligo design:* Note that the NGG PAM is **not** included in the designed oligos.



*Oligonucleotide ordering tips:* Standard de-salted oligos (usually the most inexpensive synthesis) are sufficient for cloning. If not already resuspended, dilute each oligo to 100 μM in sterile water or TE.


## Golden-Gate sgRNA cloning protocol

### 1. Oligo anneal

Component	Amount [ul]
Each oligo [100uM]	1
10X T4 ligase buffer (NEB)	1
T4 PNK (NEB)	0.5
H <sub>2</sub> O	6.5

Mix the components above and anneal in a thermal cycler with the following conditions:

37°C for 30min  
 95°C for 5 min  
 Ramp to 25°C at 5°C/min



## 2. Golden Gate reaction

Add 90ul of PCR clean H<sub>2</sub>O to the finished oligo anneal from above to dilute it 1:10. Then mix the following components:

**Note:** use **BbsI enzyme** for the non-lentiviral SAM sgRNA backbone (addgene #61424) and **BsmBI enzyme** for the lenti SAM sgRNA (zeo) backbone (addgene #61427).

Component	Amount [ul]
2X rapid ligase buffer ( <i>Enzymatics</i> )	12.5
BSA [20mg/ml] ( <i>NEB</i> )	0.125
Restriction enzyme (BbsI for <b>61424</b> or BsmBI for <b>61427</b> ) ( <i>Fermentas FD</i> )	1
T7 ligase ( <i>Enzymatics</i> )	0.125
Diluted oligo anneal (1:10) from above	1
Backbone Vector [25ng/ul]	1
H <sub>2</sub> O	9.25

Run the following program on a thermal cycler:

37C for 5 min

20C for 5 min

→ repeat for 15 cycles total

Transform 2ul of the golden gate reaction in StbI3 (or other recombination deficient) competent cells. Plate onto **Ampicillin** plates. In general, picking 2-3 colonies per guides should be sufficient to ensure a correct clone.

**Note:** it is not necessary to perform a negative control golden-gate reaction (without insert) as it will always contain colonies and is not a good indicator of cloning success.