REAL TOTAL RNA SPIN PLUS KIT

Ref. RBMER11 50 Preps.

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1. INTRODUCTION

This kit allows obtaining DNA-free Total RNA from cell cultures, animal tissues, bacteria, yeast, biologic fluids and reactions mixture using columns with a silica membrane.

For the obtaining of high quality RNA, is important to prevent the RNA degradation and to remove genomic DNA. The cells are lysed by incubation in a solution that contains chaotropic salts and immediately inactivates the RNases (which are present in all biological materials) and at the same time creates the suitable binding conditions to ease the RNA adsorption to the silica membrane. The lysis is followed by a homogenization and a filtration through a filter provided with the kit. Contaminating DNA, which is also bound to the silica membrane is removed by a DNase I solution that is directly applied onto the silica membrane. The salts, metabolites and cell components are removed by washing with two different buffers. The total RNA is eluted with nuclease-free water.

The total RNA prepared with REALTOTAL SPIN[™]PLUS[™] is useful for applications such as RT-PCR, Northern, primer extension, arrays technology and RNase protection experiments.

MUESTRA / SAMPLE	ARN PROMEDIO / RNA AVERAGE
8 x 10 ⁴ HeLa Cells	1.5 µg
4 x 10 ⁵ HeLa Cells	4 µg
1 x 10 ⁶ HeLa Cells	14µg
2 x 10 ⁶ s HeLa Cells	21 µg
2.5 x 10 ⁶ HeLa Cells	25 µg
5 x 10 ⁶ HeLa Cells	50 µg

2. COMPONENTES DEL KIT KIT COMPONENTS

	Ref.	Ref. RBMER11 50 preps	T ^a
Solución de Lisis I <i>Lysis I Solution</i>	ERC1	25 ml	RT
Solución STOP DNasa <i>DNasa STOP Solution</i>	ERC2	15 ml	RT
Solución de Lavado Concentrada <i>Concentrated Wash Solution</i>	ERC3	12.5 ml	RT
Tampón de Desalado de la Membrana <i>Membrane Desalting Buffer</i>	ERCM	25 ml	RT
Tampón de Reacción de la DNasa <i>DNase Reaction Buffer</i>	ERCD	7 ml	RT
Dnasa libre de RNasa (liofilizada) <i>RNase-free DNase (lyophilized)</i>	--	1 vial	4°C
Agua Libre de Nucleasas <i>Nuclease-free Water</i>	ERCH	15 ml	RT
Columnas de Unión ARN + Tubos Recolección <i>RNA Binding Columns + Collection Tubes</i>	ERC11	50 unid.	RT
Columnas de Filtración <i>Filtration Columns</i>	ERC12	50 unid.	RT
Tubos de Recolección <i>Collection Tubes</i>	ERC13	150 unid.	RT
Tubos de 1.5 ml <i>1.5 ml Tubes</i>	ERC14	50 unid.	RT

3. PROTOCOLO GENERAL GENERAL PROCEDURE

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3.1 Preliminary conditions

- Sample Size:**
 - Up to 5×10^6 Cells
 - Up to 30 mg tissue
 - Up to 1×10^9 bacteria
 - Up to 5×10^7 yeast
- Obtained RNA average:** Up to 70 µg
- Elution Volume:** 40 – 120 µl
- Binding Capacity:** 100 µg
- Time/ preparation:** < 30 minutes / 6 prep

It is possible to obtain up to 70 µg of total RNA from 5×10^6 cultured cells or 30 mg of tissue are used as sample. The RNA obtained can be used as template for RT-PCR. Generally, between 1-10% of the total RNA obtained from 1×10^6 cells or 10 mg of tissue is sufficient as template for RT-PCR.

The RNA prepared from such high amounts is generally DNA-free, although some DNA traces can remain in the preparation if DNA-rich materials are used. For this reason we recommend to use lower quantities of sample if it has to be used in RT-PCR.

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3.2 Preliminary Preparations

- *Reconstitute the lyophilized RNase-free DNase I, add 540 µl of nuclease-free water to the vial and incubate 1 minute at room temperature, mix until the DNase I is dissolved. Dispense into aliquots and store at -20°C. Do not thaw/freeze the aliquots more than 3 times.*
- *Add 50 ml of 100% alcohol to the concentrated wash solution, mark the date.*
- *Due to the RNases omnipresence, it is absolutely necessary that all the initial material is frozen in liquid Nitrogen and stored at -70°C. Add the lysis solution as early as possible, once it has been added the lysate can be stored at -70°C for months.*
- **ATTENTION:** *The Lysis Solution I, the STOP DNase I solution and the Membrane Desalting Buffer contains Guanidine thiocyanate. Use gloves and glasses for protection. In case of contact, wash with water.*
- *RNA is not protected against digestion until the sample material is flash frozen or disrupted in the presence of Rnase inhibiting or denaturing agents. Therefore it is important that samples are flash frozen in liquid nitrogen immediately and stored at -70°C or processed as soon possible, it is also possible to use RNAlater (AMBION). The samples can be stored in the lysis solution after the homogenization at -70°C for 1 year, at 4°C for 24 hours or up to for several hours at room temperature.*

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3.3 Protocol for RNA extraction from tissues or cultured cells

Homogenization of sample: Disrupt up to 30 mg of tissue.

IMPORTANT: *It is necessary for efficient RNA preparation that all the RNA contained in the sample is released from the cells by disruption with a mechanical homogenizator (Polytron) and that the viscosity of the sample is reduced by homogenization, be careful to keep the rotor submerged to avoid forming too much foam.*

For fresh and soft tissues use the homogenizator; for fresh and hard tissues or RNase-rich tissues grind the sample to a fine powder with liquid Nitrogen. For frozen soft tissues use the homogenizator; for the rest of tissues grind the sample to a fine powder with liquid Nitrogen.

Up to 5 x 10⁶ cultured cells can be collected by centrifugation and directly lysated by the addition of the Lysis I solution.

1. *Add 350 µl of Lysis I Solution + 3.5 µl of β-mercaptoethanol to the ground tissue or to the cell pellet and shake vigorously.*
2. *Reduce the viscosity and clear the lysate by filtration through a filtration column. Place the filtration column in a collection tube, apply the lysate and centrifuge 1 minute at 11.000 x g. Transfer the filtrate into a new centrifuge tube. Do not shake the pellet from cell debris, which can be observed at the bottom of the collection tube after the centrifugation.*
3. *Add 350 µl of 70% Ethanol to the filtrate obtained in step 2 and mix by vortex or pipette. After addition of ethanol a stringy precipitate may become visible which will not affect the RNA isolation. Be sure to load all of the precipitate on the column.*
4. *Take one RNA binding column plus collection tube and load the lysate. Centrifuge at 8.000xg (10.000 rpm.) for 30 seconds. Place the column in a new collection tube. Maximal loading capacity column is 750 µl. Repeat the procedure if larger volumes are to be processed.*
5. *Add 350 µl of Membrane Desalting Buffer and centrifuge at maximum speed (11.000 xg) for 1 minute to dry the membrane. The salts*

removal will make the following a DNase digestion much more effective. Avoid the column contact with the remaining liquid.

6. Prepare the **DNase reaction mixture** in a sterile microtube, for each isolation, add **10 µl of reconstituted DNase I to 90 µl of DNase I reaction Buffer**. Mix.
7. Apply **95 µl of the DNase reaction mixture** directly onto the centre of the silica membrane of the column. Incubate at **room temperature for 15 minutes**.
8. Add **200 µl of STOP/DNase I Solution** to the column. Centrifuge at **8.000 xg for 30 seconds**. Place the column in a new collection tube. This buffer will inactivate the DNase I.
9. Add **600 µl of Wash Solution** to the column. Centrifuge at **8.000 xg for 30 seconds**. Remove the liquid and place the column back into the collection tube.
10. Add **250 µl of Wash Solution** to the column. Centrifuge for **2 minutes at maximum speed (11.000 x g)** to completely dry the membrane. Place the column into a 1.5 ml microtube (supplied). If for any reason the liquid touches the column after the centrifugation, repeat the process.
11. Elute the RNA in **30 µl of Nuclease-free Water** (supplied), be sure that all the membrane is well soaked. **Incubate for 1 minute at room temperature**. Centrifuge at **maximum speed (11.000 x g)** for **1 minute**.
12. Add again **30 µl of Nuclease-free Water** to the same column, with the same collection tube. Centrifuge at **maximum speed (11.000 x g)** for **1 minute**. Finally it will be obtained **60 µl of RNA**

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3.4 Protocol for RNA extraction from biologicals fluids

1. It does not require the sample homogenization.
2. Add **350 µl of Lysis I Solution to 100 µl of sample**. Shake vigorously with vortex. It is not necessary to filtrate the lysate.
3. Add **350 µl of Ethanol 70 %** to the filtrate obtained in step 2 and mix with vortex or pipette. After addition of ethanol a stringy precipitate may become visible which will not affect the RNA isolation. Be sure to load all of the precipitate on the column.
4. Take one **RNA binding column** plus collection tube and load the lysate. Centrifuge at 8.000xg (**10.000 rpm.**) for **30 seconds**. Place the column in a new collection tube. Maximal loading capacity column is 750 µl. Repeat the procedure if larger volumes are to be processed.
5. Add **350 µl of Membrane Desalting Buffer** and centrifuge at **maximum speed (11.000 xg)** for **1 minute** to dry the membrane. The salts removal will make the following a DNase digestion much more effective. Avoid the column contact with the remaining liquid.
6. Prepare the **DNase reaction mixture** in a sterile microtube, for each isolation, add **10 µl of reconstituted DNase I to 90 µl of DNase I reaction Buffer**. Mix.
7. Apply **95 µl of the DNase reaction mixture** directly onto the centre of the silica membrane of the column. Incubate at **room temperature** for **15 minutes**.
8. Add **200 µl of STOP/DNase I Solution** to the column. Centrifuge at **8.000 xg** for **30 seconds**. Place the column in a new collection tube. This buffer will inactivate the DNase I.
9. Add **600 µl of Wash Solution** to the column. Centrifuge at **8.000 xg** for **30 seconds**. Remove the liquid and place the column back into the collection tube.
10. Add **250 µl of Wash Solution** to the column. Centrifuge for **2 minutes at maximum speed (11.000 x g)** to completely dry the membrane. Place the column into a 1.5 ml microtube (supplied). If for any reason the liquid touches the column after the centrifugation, repeat the process.
11. Elute the RNA in **30 µl of Nuclease-free Water** (supplied), be sure that all the membrane is well soaked. **Incubate for 1 minute at room temperature**. Centrifuge at **maximum speed (11.000 x g)** for **1 minute**.
12. Add again **30 µl of Nuclease-free Water** to the same column, with the same collection tube. Centrifuge at **maximum speed (11.000 x g)** for **1 minute**. Finally it will be obtained **60 µl of RNA**

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3.5 Protocol for RNA extraction from up to 10⁹ bacteria cells

1. *Resuspend the bacteria pellet (Gram - strains) in 100 µl TE (10 mM Tris-HCl, 1 mM EDTA; pH 8.0) containing 0.2 mg/ml lysozyme. Shake vigorously with vortex and incubate at 37°C for 10 minutes. For Gram + bacteria, resuspend the pellet in 100 µl TE (10 mM Tris-HCl, 1 mM EDTA; pH 8.0) containing 2 mg/ml lysozyme. It may be necessary to optimize the incubation time and the concentration of lysozyme depending on the bacteria strain.*
2. *Add 350 µl of Lysis I Solution + 3.5 µl of β-mercaptoethanol to the suspension and shake vigorously with vortex.*
3. *Reduce the viscosity and wash the lysate by filtration through a filtration column. Place the filtration column in a collection tube, pour the lysate and centrifuge 1 minute at 11.000 x g. Transfer the filtrate into a new centrifuge tube. Do not shake the pellet from cell debris, which can be observed at the bottom of the collection tube after the centrifugation.*
4. *Add 350 µl of Ethanol 70 % to the filtrated obtained in step 2 and mix with vortex or pipette. After addition of ethanol a stringy precipitate may become visible which will not affect the RNA isolation. Be sure to load all of the precipitate on the column.*
5. *Take one RNA binding column plus collection tube and load the lysate. Centrifuge at 8.000xg (10.000 rpm.) for 30 seconds. Place the column in a new collection tube. Maximal loading capacity column is 750 µl. Repeat the procedure if larger volumes are to be processed.*
6. *Add 350 µl of Membrane Desalting Buffer and centrifuge at maximum speed (11.000 xg) for 1 minute to dry the membrane. The salts removal will make the following a DNase digestion much more effective. Avoid the column contact with the remaining liquid.*
7. *Prepare the DNase reaction mixture in a sterile microtube, for each isolation, add 10 µl of reconstituted DNase I to 90 µl of DNase I reaction Buffer. Mix.*
9. *Apply 95 µl of the DNase reaction mixture directly onto the centre of the silica membrane of the column. Incubate at room temperature for 15 minutes.*
10. *Add 200 µl of STOP/DNase I Solution to the column. Centrifuge at 8.000 xg for 30 seconds. Place the column in a new collection tube. This buffer will inactivate the DNase I.*
11. *Add 600 µl of Wash Solution to the column. Centrifuge at 8.000 xg for 30 seconds. Remove the liquid and place the column back into the collection tube.*
12. *Add 250 µl of Wash Solution to the column. Centrifuge for 2 minutes at maximum speed (11.000 x g) to completely dry the membrane. Place the column into a 1.5 ml microtube (supplied). If for any reason the liquid touches the column after the centrifugation, repeat the process.*
13. *Elute the RNA in 30 µl of Nuclease-free Water (supplied), be sure that all the membrane is well soaked. Incubate for 1 minute at room*

temperature .Centrifuge at maximum speed (11.000 x g) for 1 minute.

14. Add again **30 µl of Nuclease-free Water** to the same column, with the same collection tube.

Centrifuge at maximum speed (11.000 x g) for 1 minute. Finally it will be obtained 60 µl of RNA

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3.6 Protocol for RNA extraction from up to 10⁸ yeast cells

1. Centrifuge **2-5 ml** of a YPD yeast culture. Resuspend the pellet **sorbitol/lyticase Buffer** (50-100 U lyticase or zymolase 1M sorbitol/ 100 mM EDTA) and incubate at **30°C for 30 minutes**. Remove the resultant spheroplastes by centrifugation (1.000 xg) for 10 minutes. It may be necessary o optimize the incubation time and the concentration of lyticase/zymolase depending on the yeast strain.
2. Add **350 µl of Lysis I Solution + 3.5 µl of β-mercaptoethanol** to the suspension and shake vigorously with vortex.
3. Reduce the viscosity and wash the lysate by filtration through a **filtration column**. Place the filtration column in a collection tube, pour the lysate and **centrifuge 1 minute at 11.000 xg**. Transfer the filtrate to a new centrifuge. **Do not shake the pellet** from cell debris which can be observed in the bottom of the collection tube after the centrifugation.

4. Add **350 µl of Ethanol 70 %** to the filtrate obtained in step 2 and mix with vortex or pipette. After addition of ethanol a stringy precipitate may become visible which will not affect the RNA isolation. Be sure to load all of the precipitate on the column.
5. Take one **RNA binding column** plus collection tube and load the lysate. Centrifuge at 8.000xg (**10.000 rpm.**) for **30 seconds**. Place the column in a new collection tube. Maximal loading capacity column is 750 µl. Repeat the

procedure if larger volumes are to be processed.

6. Add **350 µl of Membrane Desalting Buffer** and centrifuge at **maximum speed (11.000 xg)** for **1 minute** to dry the membrane. The salts removal will make the following a DNase digestion much more effective. Avoid the column contact with the remaining liquid.
7. Prepare the **DNase reaction mixture** in a sterile microtube, for each isolation, add **10 µl of reconstituted DNase I** to **90 µl of DNase I reaction Buffer**. Mix.
8. Apply **95 µl of the DNase reaction mixture** directly onto the centre of the silica membrane of the column. Incubate at **room temperature** for **15 minutes**.
9. Add **200 µl of STOP/DNase I Solution** to the column. Centrifuge at **8.000 xg** for **30 seconds**. Place the column in a new collection tube. This buffer will inactivate the DNase I.

10. Add **600 µl of Wash Solution** to the column. Centrifuge at **8.000 xg** for **30 seconds**. Remove

the liquid and place the column back into the collection tube.

11. Add **250 µl of Wash Solution** to the column. Centrifuge for **2 minutes at maximum speed (11.000 x g)** to completely dry the membrane. Place the column into a 1.5 ml microtube (supplied). If for any reason the liquid touches the column after the centrifugation, repeat the process.
12. Elute the RNA in **30 µl of Nuclease-free Water** (supplied), be sure that all the membrane is well soaked. **Incubate for 1 minute at room temperature**. Centrifuge at **maximum speed (11.000 x g)** for **1 minute**.
13. Add again **30 µl of Nuclease-free Water** to the same column, with the same collection tube. Centrifuge at **maximum speed (11.000 x g)** for **1 minute**. Finally it will be obtained **60 µl of RNA**

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3.7 Protocol for RNA Clean-up from reaction mixtures

1. *The homogenization is not necessary.*
2. *Add 2.5 volumes of Lysis I Solution per 1 volume of sample and shake with vortex. . Maximal loading capacity column is 750 µl. Repeat the procedure if larger volumes are to be processed.*
3. *The filtration is not necessary.*
4. *Add the same quantity of 70% ethanol to the lysate and shake with vortex.*
5. *Take one RNA binding column plus collection tube and load the lysate. Centrifuge at 8.000xg (10.000 rpm.) for 30 seconds. Place the column in a new collection tube. Maximal loading capacity column is 750 µl. Repeat the procedure if larger volumes are to be processed.*
6. *Add 350 µl of Membrane Desalting Buffer and centrifuge at maximum speed (11.000 xg) for 1 minute to dry the membrane. The salts removal will make the following a DNase digestion much more effective. Avoid the column contact with the remaining liquid.*
7. *Prepare the DNase reaction mixture in a sterile microtube, for each isolation, add 10 µl of reconstituted DNase I to 90 µl of DNase I reaction Buffer. Mix.*
8. *Apply 95 µl of the DNase reaction mixture directly onto the centre of the silica membrane of the column. Incubate at room temperature for 15 minutes.*
9. *Add 200 µl of STOP/DNase I Solution to the column. Centrifuge at 8.000 xg for 30 seconds. Place the column in a new collection tube. This buffer will inactivate the DNase I.*
10. *Add 600 µl of Wash Solution to the column. Centrifuge at 8.000 xg for 30 seconds. Remove the liquid and place the column back into the collection tube.*
11. *Add 250 µl of Wash Solution to the column. Centrifuge for 2 minutes at maximum speed (11.000 x g) to completely dry the membrane. Place the column into a 1.5 ml microtube (supplied). If for any reason the liquid touches the column after the centrifugation, repeat the process.*
12. *Elute the RNA in 30 µl of Nuclease-free Water (supplied), be sure that all the membrane is well soaked. Incubate for 1 minute at room temperature .Centrifuge at maximum speed (11.000 x g) for 1 minute.*
13. *Add again 30 µl of Nuclease-free Water to the same column, with the same collection tube. Centrifuge at maximum speed (11.000 x g) for 1 minute. Finally it will be obtained 60 µl of RNA*

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1. Degraded RNA or not RNA obtaining

1.1. *Try to create an RNase-free area inside your lab. You can use RNase Remove (Ref.RD056) for cleaning the material and surfaces. Change gloves frequently. It is recommended to use sterile tubes and to keep them close when possible during the process.*

2. Low quality or low quantity of RNA obtained

2.1. *Reagents not correctly reconstituted. Add the indicated quantity of nuclease-free water to the DNase I vial and Ethanol to the concentrated wash Solution.*

2.2. *The samples and the reagents have not been completely mixed.*

2.3. *The ethanol has not been added after the lysis. The RNA will only bind to the silica membrane under the presence of alcohol*

2.4. *Store the kit components as indicated in the protocol. Low temperature storages can cause salts precipitation.*

2.5. *Keep the containers correctly closed to prevent evaporation*

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- 2.6. *The sample hasn't been correctly preserved. Use, if possible, fresh material. If it is not possible, freeze the sample in liquid Nitrogen and store them at -70°C until you add the Lysis Solution*
- 2.7. *The initial material hasn't been enough homogenized. Always spray the samples in presence of liquid Nitrogen. Use the filtration columns.*
3. **Overloaded membrane**
 - 3.1. *Too much initial material has been used. An overloaded membrane can lead to a low yield. Reduce the sample quantity or use a bigger quantity of Lysis Solution.*
4. **Contamination with Genomic DNA**
 - 4.1. *Not active DNase I. Reconstitute and store the DNase I according to the supplied instructions.*
 - 4.2. *Pipette the DNase I solution directly on the centre of the silica membrane.*
 - 4.3. *Reduce the sample starting quantity.*
5. **Low RNA yield in following applications or experiments**
 - 5.1. *The Wash solution that contains ethanol has not been completely removed. Be sure that the end of the column does not touch the liquid after 2nd wash.*
 - 5.2. *Check that the Wash solution has been balanced at room temperature before use. Washing at low temperatures produces a decrease on the salt removal efficiency.*

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BIBLIOGRAPHY

- Martin, A., Adam, H., Díaz-Mendoza, M., Żurczak, M., González-Schain, N. D., & Suárez-López, P. (2009). Graft-transmissible induction of potato tuberization by the microRNA miR172. *Development*, 136(17), 2873-2881.
- Escames, G., López, L. C., Ortiz, F., López, A., García, J. A., Ros, E., & Acuña-Castroviejo, D. (2007). Attenuation of cardiac mitochondrial dysfunction by melatonin in septic mice. *FEBS Journal*, 274(8), 2135-2147.