

# LC3B Antibody Kit for Autophagy

Catalog no. L10382

**Table 1.** Contents and storage information.


Material	Amount	Concentration	Storage	Stability
LC3B, rabbit polyclonal antibody – unconjugated (Component A)	55 µL	1 mg/mL in 10 mM PBS, pH 7.4, 15 mM azide	<ul style="list-style-type: none"> <li>• ≤-20°C</li> <li>• Protect from light</li> </ul>	When stored as directed the product is stable for at least 6 months from the date of receipt.
Chloroquine diphosphate (Component B)	1 mL	30 mM aqueous solution		

## Introduction


Autophagy describes the segregation and delivery of cytoplasmic cargo, including proteins and organelles, for degradation by hydrolytic enzymes in “autophagolysosomes”, also referred to as “autolysosomes”. Although first described in 1963, it has only been in the past decade that this pathway has been the subject of intense research to gain further insight into the role basal autophagy plays in cell homeostasis and development. Efforts are also directed to further elucidate the role of induced autophagy as a cell survival response to stress, microbial infection, and disease (e.g., neurodegeneration, cancer).<sup>1-3</sup>

The LC3B protein plays a critical role in autophagy. Normally, this protein resides in the cytosol, but following cleavage and lipidation with phosphatidylethanolamine, LC3 associates with the phagophore. This localization can be used as a general marker for autophagic membranes (Figure 1).

Each LC3B Antibody Kit for Autophagy includes a rabbit polyclonal antibody against LC3B that has been validated for use in fluorescence microscopy and high content imaging and analysis. The kit also includes chloroquine diphosphate for artificially generating autophagosomes (Figure 2). Following chloroquine diphosphate treatment, lysosomal pH increases and the normal autophagic flux is disrupted, resulting in autophagosome accumulation.



**Figure 1.** Schematic depiction of the autophagy pathway in a eukaryotic cell. The first step involves the formation and elongation of the isolation membranes or phagophore. The second step entails the expansion and sequestering of the cytoplasm and formation of the double-membrane autophagosome and includes the association of the cytosolic LC3 protein. Fusion of lysosomes with autophagosome to generate the autolysosome is the penultimate step. In the fourth and final phase, the cargo is degraded.



**Figure 2.** Autophagy detection with the LC3B Antibody Kit. HeLa, A549 and HepG2 cells were treated with 50  $\mu$ M chloroquine for 16 hours at 37°C. Following fixation and permeabilization, autophagosomes were stained with LC3B rabbit polyclonal antibody and visualized with Alexa Fluor® 647 goat anti-rabbit IgG. The fluorescence intensity of the autophagosomes and the cytosol were quantified using Slidebook™ digital microscopy software.

## Before Starting

---

### Materials Required but Not Provided

- Phosphate buffered saline (PBS), pH 7.2–7.6
- Fixative (i.e., 3.7% formaldehyde in PBS)
- Permeabilization reagent (i.e., 0.2% Triton® X-100 in PBS)
- Blocking buffer
- Anti-rabbit secondary detection reagent

## Experimental Protocols

---

The following protocol was optimized for fluorescence microscopy and high-content imaging and analysis using cells fixed with 3.7% formaldehyde in PBS and permeabilized with 0.2% Triton® X-100 in PBS, but the LC3B antibody is also amenable to other fixation/permeabilization reagents.

- 1.1** *Optional:* Treat cells with 30–100  $\mu$ M chloroquine (Component B) for 12–16 hours.
- 1.2** Dilute the LC3B rabbit polyclonal antibody in blocking buffer to prepare 0.5  $\mu$ g/mL working solution. Note that this final concentration may require further optimization.
- 1.3** Add 3.7% formaldehyde in PBS to the cells and incubate for 15 minutes at room temperature.
- 1.4** Remove the fixative and wash the cells three times with PBS.
- 1.5** Remove the wash solution, add 0.2% Triton® X-100 in PBS to the cells, and incubate them for 15 minutes at room temperature.
- 1.6** Remove the permeabilization buffer, add the diluted primary antibody (prepared in Step 1.2) to the cells, and incubate them for 1 hour at room temperature.
- 1.7** Remove the diluted primary antibody and wash the cells three times with PBS.
- 1.8** Remove the wash solution and incubate the cells with an anti-rabbit secondary antibody.
- 1.9** Wash the cells with PBS and perform any additional staining (i.e., DNA counterstain, other antibodies).
- 1.10** Image the cells using the appropriate filters for the anti-rabbit secondary detection reagent. Autophagosomes are vesicular structures typically located in the perinuclear region.

## References

---

- 1.** Molecular Cell Biology 8, 931 (2007); **2.** Genes & Development 21, 2861 (2007); **3.** Drug Discovery 6, 304 (2007).

**Product List** Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
L10382	LC3B Antibody Kit for Autophagy *rabbit polyclonal LC3B* *includes autophagosome inducer*	1 kit
<b>Related Products</b>		
P36235	Premo™ Autophagy Sensor LC3B-GFP	1 kit
P36236	Premo™ Autophagy Sensor LC3B-RFP	1 kit
A11008	Alexa Fluor® 488 goat anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A11012	Alexa Fluor® 594 goat anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A11034	Alexa Fluor® 488 goat anti-rabbit IgG (H+L) *highly cross-adsorbed* *2 mg/mL*	0.5 mL
A11037	Alexa Fluor® 594 goat anti-rabbit IgG (H+L) *highly cross-adsorbed* *2 mg/mL*	0.5 mL
A11039	Alexa Fluor® 594 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A11046	Alexa Fluor® 350 goat anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21068	Alexa Fluor® 350 goat anti-rabbit IgG (H+L) *highly cross-adsorbed* *2 mg/mL*	0.5 mL
A21206	Alexa Fluor® 488 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21207	Alexa Fluor® 594 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21244	Alexa Fluor® 647 goat anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21245	Alexa Fluor® 647 goat anti-rabbit IgG (H+L) *highly cross-adsorbed* *2 mg/mL*	0.5 mL
A21248	Alexa Fluor® 555 goat anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21249	Alexa Fluor® 555 goat anti-rabbit IgG (H+L) *highly cross-adsorbed* *2 mg/mL*	0.5 mL
A21441	Alexa Fluor® 488 chicken anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21442	Alexa Fluor® 594 chicken anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21443	Alexa Fluor® 647 chicken anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A31572	Alexa Fluor® 555 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A31573	Alexa Fluor® 647 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL

**Contact Information**

**Molecular Probes, Inc.**

29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

**Customer Service:**

6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probesorder@invitrogen.com

**Toll-Free Ordering for USA:**

Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

**Technical Service:**

8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probestech@invitrogen.com

**Invitrogen European Headquarters**

Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

For country-specific contact information, visit [www.invitrogen.com](http://www.invitrogen.com).

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

**Limited Use Label License No. 223: Labeling and Detection Technology**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2010, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.