

Novus Biologicals uses cookies to provide you with a great website experience. By continuing to use this website you acknowledge this and agree to our [cookie policy](#). [OK](#) [Learn More](#)

Support: 1-813-300-0220 | info@novusbio.com

[Items in Cart \(0\)](#) [Quick Order](#) [Newsletter](#)

[sign in or register](#)



[Products](#)

[Resources](#)

[Research Areas](#)

[Customer Care](#)

[About Us](#)

[Contact](#)

[Home](#) » [Support](#) » Protocol: Sandwich ELISA (Colorimetric) – Direct Detection

Protocol: Sandwich ELISA (Colorimetric) – Direct Detection

Related Links

[ELISA Protocols, Troubleshooting, & Scientific Resources](#)

[ELISA Metrics Infographic](#)

[View all protocols](#)

Capture and detection antibodies should be matched for an ELISA. They should not bind the same epitope or recognize epitopes in close proximity.

MATERIALS

- **1X Phosphate Buffered Saline (PBS):** Dissolve 1.44 g Na_2HPO_4 (10 mM), 0.24 g KH_2PO_4 (1.8 mM), 8 g NaCl (137 mM), and 0.2 g KCl (2.7 mM) in 800mL of dH_2O . Adjust the pH to 7.4 with 0.1N HCl, and add dH_2O up to 1L.
- **Coating Buffer:** Dissolve anhydrous 5.3 g Na_2CO_3 (50mM) and anhydrous 4.2 g NaHCO_3 , (50mM) in dH_2O . Adjust pH to 9.6 and add dH_2O up to 1L.
- **Reagent Diluent:** Dissolve BSA in Phosphate Buffered Saline (1X PBS) to 1% concentration (w/v).
- **Wash buffer:** 1X PBS containing 0.05% v/v Tween-20
- **Capture Antibody Solution:** Add the primary capture antibody to the Coating Buffer. If a working antibody concentration is not available, use the below table for recommended starting concentrations. The clonality of the antibody will impact the suitable concentration range.

	Monoclonal Capture/ Polyclonal Detection	Monoclonal Capture/ Monoclonal Detection	Polyclonal Capture/ Polyclonal Detection
Capture Antibody	1-8 $\mu\text{g}/\text{mL}$	0.5-4 $\mu\text{g}/\text{mL}$	0.2-0.8 $\mu\text{g}/\text{mL}$
Detection Antibody	50-400 ng/mL	0.25-2 $\mu\text{g}/\text{ml}$	50-400 ng/mL

- **Detection Antibody Solution:** Add the enzyme conjugated detection antibody to the Reagent Diluent. Common enzyme conjugates for ELISAs include Horse Radish Peroxidase (HRP) and Alkaline Phosphatase (AP). If a working antibody concentration is not available, use the above table for recommended starting concentrations. The clonality of the antibody will impact the suitable concentration range.
- **Substrates and Stop Solutions:** Follow datasheets for preparation of Substrate Solutions.

Conjugates	Substrates	Stop Solution	Absorbance Settings
HRP	TMB (3,3',5,5'-Tetramethylbenzidine)	2M H_2SO_4	450nm
AP	PNPP (Para-Nitro Phenyl Phosphate)	1M NaOH	405nm

SAMPLE PREPARATION

Dilute the experimental samples (serum, plasma, lysates or other biofluids) and perform a serial dilution of the standard (purified form of the analyte under consideration) in Reagent Diluent. If available, use the literature to determine the optimal concentration and linear working range. Samples should be prepared in duplicate, or preferably in triplicate.

METHODS

1. Use an ELISA compatible PVC clear microtiter plate (96-well) and determine the number of wells needed for the assay.
2. Dispense 100 μl of Capture Antibody Solution into the wells. Apply sealing tape to the top of the plate to prevent evaporation. Incubate the plate overnight at 4°C. After incubation remove the

Novus Biologicals uses cookies to provide you with a great website experience. By continuing to use this website you acknowledge this and agree to our [cookie policy](#). [OK](#) [Learn More](#)

out the solution. Wash the plate 2 more times and pat the inverted plate on a paper towel to dry.

4. Block plates by dispensing 300 μ l of Reagent Diluent into each well and incubate the plate for 60 minutes at room temperature. After incubation, remove the tape. Invert the plate and flick out the solution.
5. Repeat step 3.
6. Add 100 μ l of diluted samples and standards to the appropriate wells and cover the plate.
7. Incubate the plate for 2 hours at room temperature. After incubation, uncover the plate and aspirate each well.
8. Repeat step 3.
9. Dispense 100 μ l of Detection Antibody Solution to each well. Incubate for 2 hours at room temperature and cover the plate. After incubation, uncover the plate. Invert the plate and flick out the solution.
10. Repeat step 3.
11. Dispense 100 μ l of the Substrate Solution into each well. Cover the plate and incubate for 20 to 30 minutes at room temperature. Avoid placing the plate in direct sunlight.

Tip: If no color develops within 30 minutes, incubation time can be extended until an appropriate color change/signal is acquired. TMB substrate reacts with HRP to generate a blue colored product which changes to yellow upon the addition of Stop Solution, whereas the PNPP substrate reacts with ALP to generate a yellow colored product.

12. After color development, remove the cover and dispense 50 μ l of Stop Solution into each well to stop the enzymatic reaction.
13. Immediately measure the absorbance of each well using a spectrophotometer/plate reader with the appropriate absorbance setting.

Tip: The presence of air bubbles in wells may impact the results and it is advisable to spin the plate in a centrifuge.

Popular Products

[Antibodies](#)
[Secondary Antibodies](#)
[Isotype Controls](#)
[Proteins and Peptides](#)
[ELISA Kits](#)
[Sample Size](#)
[New Products](#)

Company Information

[About Us](#)
[Careers](#)
[Novus Guarantee](#)
[Contact Us](#)
[Our Team](#)

Support

[Technical Support](#)
[FAQs](#)
[BioInformatics](#)
[Help Finding Antibodies](#)
[Protocols](#)

Stay Connected

[Newsletter Signup](#)

[Find us on Facebook](#)

[Follow us on Twitter](#)

[Watch us on YouTube](#)

[Connect with us on LinkedIn](#)

[Connect with us on Pinterest](#)

[Connect with us on Instagram](#)

Europe Office Bio-Techne Ltd. Abingdon Oxon
Phone: +31 800 0225607 | info.emea@bio-techne.com | Fax: +31 2 071 68337
Office Hours: 8am-5pm GMT

[Visit our other sites](#)