

abcam

# Native chromatin immunoprecipitation protocol

Detailed procedure and tips for ChIP

## The preparation of native chromatin from cultured human cells.

All solutions need to be ice cold.

Sucrose containing solutions must be made up fresh on the day.

**Add protease inhibitors to all lysis solutions before use (0.1mM PMSAF and complete mini protease inhibitors; commercially available).**

1. Cultured cells (e.g. HL-60 or lymphoblastoids) are grown to a density of approximately  $1 \times 10^6$  cells/ml, until they are in log phase.
2. Harvest cells: centrifuge samples (1,000 g, 10 min, 4°C) and wash the cell pellet 3 x with ice cold PBS (Phosphate buffered saline).

**It is essential that 5 mM Na butyrate is present in all solutions throughout chromatin isolation when using antibodies to acetylated histones to prevent deacetylation.**

3. Resuspend cell pellet in TBS (Tris buffered saline) at  $2 \times 10^7$  cells/ml and add an equal volume of 1.0% v/v Tween 40 in TBS. Add PMSF to a final concentration of 0.5 mM. Leave stirring gently on ice for 1 hr (Transfer the suspension into a 50 ml tube with a small magnetic bar or flea; place the tube in ice on top of a magnetic stirrer).
4. Transfer cell lysate to an all-glass homogenizer and homogenize 7 ml aliquots with seven strokes using an 'A' or 'tight' pestle. Check that nuclei have been released by phase-contrast microscopy; intact cells should have the central dark region of the nucleus surrounded by a halo, which is the less dense cytoplasm.

**You may have to increase or decrease this homogenization step to maximize the yield of nuclei depending on cell line. Keep cells on ice between the rounds of homogenisations.**

5. Centrifuge samples (8,000 g, 20 min, 4°C).
6. Resuspend nuclei pellet in 25% [w/v] sucrose/TBS at  $4 \times 10^6$  nuclei/ml and underlay with 0.5 vol of 50% [w/v] sucrose/ TBS; centrifuge the samples (10,000 g, 15 min, 4°C).
7. Discard supernatant and wash nuclei pellet in 5 ml 25% [w/v] sucrose/TBS; centrifuge samples (10,000 x g, 15 min, 4°C).
8. Resuspend nuclei pellet in 5 ml digestion buffer and check absorbance ratios at 260 nm and 280 nm for a diluted sample of the nuclei suspension; calculate the approximate DNA concentration from the A260 reading (the ratio of A260/A280 should be about 1.1).

**The yield of chromatin (in  $\mu\text{g}$ ) is given by:  $A_{260} \times \text{dilution factor} \times \text{volume} \times 50$ .  
Centrifuge samples (10,000 rpm, 10 min, 4°C) and resuspend the nuclei pellet at 0.5 mg/ml in 1.7 ml Eppendorf tube(s). Divide into 1 ml aliquots, if necessary.**

## Micrococcal nuclease digestion

Normally we add 50 U micrococcal nuclease per 0.5 mg DNA, in a reaction volume of 1.0 ml. This is usually provided as a powder; dissolve the micrococcal nuclease in dH<sub>2</sub>O to the required concentration and store as small aliquots at -20°C. Aliquots may be re-frozen and re-used once. This step needs to be carefully controlled, especially in the initial preparations.

**High concentrations of micrococcal nuclease may over-digest the chromatin, leading to sub-nucleosomal particles. You should aim to obtain a long/medium oligonucleosome ladder. If pure mononucleosome preparations are required carry out a linear sucrose gradient (5-20%), this will increase resolution.**

1. Perform micrococcal nuclease digestions at 37°C for 5 min.
2. Stop reaction by addition of 0.2 M EDTA to a final concentration of 5 mM.
3. Place all samples on ice for 5 min; centrifuge samples (12,000 g, 5 min).
4. Remove and keep the first S/N (this is called the S1 fraction; total vol 1.0 ml); store overnight at 4°C.
5. Resuspend the pellet in 1.0 ml Lysis buffer and dialyze overnight against 2 l of the same buffer.
6. After overnight dialysis centrifuge samples (1,800 g, 10 min, 4°C).
7. Remove and keep the supernatant (called the S2 fraction; total vol about 1.2 ml after dialysis); store at 4°C.
8. Resuspend insoluble pelleted material in 200 µl lysis buffer (called the P fraction).

## Analysis of soluble chromatin fractions

1. Check A<sub>260</sub>/A<sub>280</sub> in all samples; the ratios for S1, S2 and P fractions are approximately 1.7, 1.5 and 1.3 respectively. Analyze all samples by 1.2% agarose gel electrophoresis.

**Do not place ethidium bromide in the agarose gel or the electrophoresis buffer, because of the presence of SDS (see below).**

Preparation of samples: x µl (total of 5 µg) chromatin fraction (S1, S2 and P), y µl dH<sub>2</sub>O (x+y = 25 µl), 3 µl 1% [w/v] SDS (final conc 0.1%), 2 µl gel loading buffer, containing bromophenol blue.

2. Stain the gel with 0.5 µg/ml ethidium bromide after the run has finished.

## Immunoprecipitation

1. 100-200 µg unfixed chromatin + 100-200 µl affinity purified antibody (50-100 µg Ig) and the final volume made up to 1.0 ml with incubation buffer. A negative control, with no added antibody, also needs to be set up to test for any nonspecific binding of the chromatin to the protein A Sepharose.
2. Incubate overnight at 4°C on a slow rotating turntable. Add 200 µl 50% v/v protein A Sepharose; use a siliconized pipette with the tip cut off to make this step easier. Incubate for 3 hr at room temperature on a fast rotating turntable. (Make sure that the Sepharose is in a suspension at all times).
3. Centrifuge samples (2,000 g, 10 min, 4°C), remove and keep the S/N; this is the unbound (or "U") fraction.
4. Resuspend the Sepharose pellet in 1 ml buffer A and layer onto 9 ml of the same buffer using a siliconized pasteur pipette and siliconized 15 ml tube.
5. Centrifuge samples (2,000 g, 10 min, 4°C), discard the S/N and wash the Sepharose sequentially in 10 ml buffer B and buffer C.
6. Finally, resuspend the Sepharose in 1 ml buffer C and transfer back to siliconized Eppendorfs.
7. Centrifuge samples (2,000 g, 10 min, 4°C) and resuspend the sepharose pellet in 250 µl 1.0% SDS / incubation buffer and incubate for 15 min at RT on a fast turntable. (Ensure that the Sepharose is thoroughly resuspended at all times).
8. Centrifuge the samples (2,000 g, 10 min, 4°C) and remove and keep S/N; this is the bound (or "B") fraction.
9. Wash the Sepharose in 250 µl 1.0% SDS / incubation buffer and centrifuge immediately (2,000 g, 10 min, 4°C). Remove the S/N and pool with the previous bound fraction from the previous step.

## DNA Isolation

**Add 500 µl incubation buffer to each bound fraction, to reduce the SDS concentration to 0.5% SDS. Unbound and bound fractions then treated as follows:**

1. Add 0.33 vol (330 µl) phenol/chloroform; vortex and spin (13,000 rpm, 10 min, microcentrifuge). Keep the organic phase and interface; this is used to isolate immunoprecipitated proteins (see below).
2. Transfer the aqueous supernatant to an equal volume (1.0 ml) of phenol/chloroform; vortex and spin (13,000 g, 10 min, microcentrifuge).
3. Transfer supernatant to an equal volume (1.0 ml) of chloroform; vortex and spin (13,000 g, 10 min, microcentrifuge).
4. Transfer S/N to a clean centrifuge tube and add 0.1 vol (100 µl) 4 M LiCl, 50 µg glycogen (Molecular biology grade, dissolved in dH<sub>2</sub>O at 2 mg/ml) as a carrier and 4 vol of ethanol. Vortex thoroughly and leave at -20°C overnight.
5. Centrifuge samples (13,000 g, 15 min) to precipitate the DNA.
6. Wash the pellet with 70% ethanol (molecular biology grade) and redissolve the DNA in 250 µl TE buffer.
7. Store samples at -20°C or proceed with detection method (PCR, microarray, etc).

8. PCR is used to quantify DNA levels of specific loci. This is analyzed semi-quantitatively (analyses of PCR end-product by agarose gel) using primers which can be designed using this tool.

**Alternatively, DNA levels are quantitatively measured by real-time PCR. Primers and probes are often designed using software provided with the real-time PCR apparatus.**

## Protein Isolation

1. To the first phenol/chloroform phase (see DNA isolation; step 1) add 5  $\mu$ l of a 1 mg/ml solution of BSA (to be used as a carrier), 0.01 vol (4  $\mu$ l) 10 M H<sub>2</sub>SO<sub>4</sub> and 12 vol of acetone.
2. After overnight precipitation at -20°C, wash the protein pellets once in acidified acetone (1:6 100 mM H<sub>2</sub>SO<sub>4</sub>:acetone) and 3 times in dry acetone. Proteins can be analyzed by SDS-PAGE.

## Solutions

### 10 x TBS

0.1 M Tris-HCl (pH 7.5)  
1.5 M NaCl  
30 mM CaCl<sub>2</sub>  
20 mM MgCl<sub>2</sub>  
50 mM Na butyrate (pH 8.0)

### Digestion buffer

0.32 M sucrose  
50 mM Tris-HCl (pH 7.5)  
4 mM MgCl<sub>2</sub>  
1 mM CaCl<sub>2</sub>  
0.1 mM PMSF  
5 mM Na butyrate

### Lysis buffer

1.0 mM Tris-HCl (pH 7.4)  
0.2 mM Na<sub>2</sub>EDTA  
0.2 mM PMSF  
5 mM Na butyrate

### Incubation buffer

50 mM NaCl  
20 mM Tris-HCl (pH 7.5)  
20 mM Na butyrate  
5 mM Na<sub>2</sub>EDTA  
0.1 mM PMSF

### Buffer A

50 mM Tris-HCl, (pH 7.5)  
10 mM EDTA  
5 mM Na butyrate  
50 mM NaCl

**Buffer B**

50 mM Tris-HCL (pH 7.5)  
10 mM EDTA  
5 mM Na butyrate  
100 mM NaCl

**Buffer C**

50 mM Tris-HCL (pH 7.5)  
10 mM EDTA  
5 mM Na butyrate  
150 mM NaCl

**Protein A Sepharose**

Pre-swell protein A Sepharose overnight in buffer A at 4°C. Centrifuge (10,000 x g, 10 min) and resuspend pellet in approximately an equal volume (50% v/v) of buffer A.

*(Adapted from protocols used by Laura O'Neill and Prof. Bryan Turner. University of Birmingham)*