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INTRODUCTION

The peqGOLD family of products is an innovative system that radically simplifies extraction and purification of nucleic acids from a variety of sources. Key to the system is the PerfectBind matrix that avidly, but reversibly, binds DNA or RNA under certain optimal conditions allowing proteins and other contaminants to be removed. Nucleic acids are easily eluted with deionized water or low salt buffer.

PerfectBind DNA Columns facilitate the binding, washing and elution steps thus enabling multiple samples to be simultaneously processed. Yields vary according to plasmid copy number, *E.coli* strain and conditions of growth, but 5 ml of overnight culture with high copy-number plasmids in LB medium typically produce 15 to 25 µg plasmid DNA. For more DNA, we recommend the peqGOLD Plasmid Miniprep Kit II that yields 40 to 75 µg DNA from 10 to 15 ml culture when using high copy-number plasmids.

Plasmid-DNA isolated with the peqGOLD Plasmid Miniprep Kit I is suitable for automated fluorescent DNA sequencing, restriction endonuclease digestion and other applications.

peqGOLD Plasmid Mini Kits I are available with Safety-Line or Classic-Line columns (Safety-Line, Ordering No. 12-6943-xx; Classic-Line, Ordering No. 12-6942-xx). Safety-Line columns can be closed tightly by lids to avoid cross-contamination more effectively. Classic-Line columns do not have lids for a more comfortable handling.

PRINCIPLE

The Plasmid Miniprep Kits combine the power of PerfectBind technology with the time-tested consistency of alkaline-SDS lysis of bacterial cells to deliver high quality DNA. Plasmid DNA binds effectively to the silica membrane of the PerfectBind DNA Column and can easily be purified from contaminants and enzyme inhibitors. Purified plasmid DNA is simply eluted from the column with Elution Buffer or water.

KIT COMPONENTS

peqGOLD Plasmid Miniprep Kit I	5 Purifications	50 Purifications	200 Purifications
Order No. Classic-Line	12-6942-00	12-6942-01	12-6942-02
Order No. Safety-Line	12-6943-00	12-6943-01	12-6943-02
Components			
PerfectBind DNA Columns	5	50	200
2.0 ml Collection Tubes	5	50	200
Solution I	5 ml	20 ml	60 ml
Solution II	5 ml	20 ml	60 ml
Solution III	5 ml	20 ml	80 ml
PW Plasmid Buffer	5 ml	30 ml	110 ml
DNA Wash Buffer (conc.)	12 ml	40 ml	3 x 40 ml
RNase A	25 µl	70 µl	200 µl
Elution Buffer (5 mM Tris-HCl, pH 8.5)	1.5 ml	15 ml	60 ml
Instruction manual	1	1	1

STORAGE AND STABILITY

All peqGOLD Plasmid Miniprep Kit I components are stable for at least 12 months from the date of purchase. RNase A and solution I completed with RNase should be stored at 4 °C, all other components should be stored at room temperature (22 °C – 25 °C).

BEFORE STARTING

Briefly examine this booklet and become familiar with each step. Prepare all components and have the necessary materials ready before starting.

- ! Before using, Solution I must be completed with RNase A. Add the complete volume of RNase A to the bottle of Solution I and mix carefully. Store the completed Solution I at 4 °C.
- ! Solution II should be closed firmly when not used.
- ! DNA Wash Buffer is concentrated and has to be completed with absolute ethanol as follows:

12-6942-00	Mix 12 ml DNA Wash Buffer with 18 ml 100 % EtOH.
12-6943-00	
12-6942-01	Mix 40 ml DNA Wash Buffer with 60 ml 100 % EtOH.
12-6943-01	
12-6942-02	Mix 3 x 40 ml DNA Wash Buffer with 3 x 60 ml 100 % EtOH.
12-6943-02	

- ! Store diluted DNA Wash Buffer at room temperature.
- ! All steps must be carried out at room temperature.

DANGEROUS COMPONENTS

peqGOLD Plasmid Miniprep I Kit

Components	Signal word / symbols	Dangerous components	H and P statements
PerfectBind DNA Columns	-	-	-
2.0 ml Collection Tubes	-	-	-
Solution I	-	-	-
Solution II	DANGER 	sodium dodecyl sulphate 1-≤2.5%, sodium hydroxide 0.5- <1%, CAS: 151-21-3, 1310-73-2	H314, P101, P102, P103, P260, P280, P303+P361+P353, P305+P351+P338, P405, P501
Solution III	DANGER  	guanidinium chloride 25-50%, acetic acid 10- <25%, CAS: 50-01-1, 64-19-7	H302, H315, H318, P101, P102, P103, P280, P264, P305+P351+P338, P362, P403+P233, P501
PW Plasmid Buffer	DANGER  	guanidinium chloride 25-50%, propan-2-ol 25-50%, CAS: 50-01-1, 67-63-0	H225, H302, H315, H319, H336, P101, P102, P103, P210, P261, P303+P361+P353, P305+P351+P338, P405, P501
DNA Wash Buffer (Konz.)	-	-	-
RNase A	DANGER  	Nuclease, ribo- 20 mg/ml, CAS: 9001-99-4	H317, H334, P261S, P280sh, P302+P352, P304+P340, P333+P313, P342+P311, P363
Elution Buffer (5 mM Tris-HCl, pH 8.5)	-	-	-

H and P statements	Descriptions
H225	Highly flammable liquid and vapour.
H302	Harmful if swallowed.
H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H336	May cause drowsiness or dizziness.
P101	If medical advice is needed, have product container or label at hand.
P102	Keep out of reach of children.
P103	Read label before use.
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P260	Do not breathe dust/fume/gas/mist/vapours/spray.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P261S	Avoid breathing dust.
P264	Wash ... thoroughly after handling.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P280sh	Wear protective gloves/eye protection.
P302+P352	IF ON SKIN: Wash with plenty of water/...
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P342+P311	If experiencing respiratory symptoms: Call a POISON CENTER/doctor/...
P362	Take off contaminated clothing.
P363	Wash contaminated clothing before reuse.
P403+P233	Store in a well-ventilated place. Keep container tightly closed.
P405	Store locked up.
P501	Dispose of contents/container to ...

PEQGOLD PLASMID MINIPREP KIT I PROTOCOL

Materials required, to be supplied by the user:

- ! Appropriate Bacterial Growth Medium
- ! 100 % Ethanol
- ! Sterile deionized water (optional)
- ! Sterile pipet tips and centrifuge tubes

A. High Copy-Number Plasmids

1. Bacterial Culture

Inoculate 5 ml growth medium containing the appropriate antibiotic in a 10 - 20 ml culture tube with *E.coli* carrying the desired plasmid and grow at 37 °C with shaking for 12 to 16 h.

Pellet overnight culture by centrifugation for 1 min at 10.000 x g* in a 1.5 ml centrifuge tube or 10 min at 5.000 x g* in a 15 ml centrifuge tube. Pour off supernatant and discard.

For High-Copy Number Plasmids, do not use more than 5 ml overnight culture, otherwise the column may become overloaded. If bigger volumes have to be processed, we recommend processing of multiple samples from the same culture. Alternatively, work with the peqGOLD Plasmid Miniprep Kit II (product no. 12-6946 / 12-6945), that allows processing of 10 - 15 ml culture using the mini-column format and generally yields 40 - 75 µg plasmid DNA.

2. Lysis of Bacteria

Add 250 µl Solution I/RNase A to the bacterial pellet and completely resuspend cells by vortexing.

Complete resuspension of cell pellet is essential for obtaining high yields.

Add 250 µl Solution II and gently mix by inverting and rotating tube 6 - 10 times to obtain a cleared lysate. A 2 min incubation at room temperature may be necessary.

Avoid vigorous mixing as this will shear chromosomal DNA and lower plasmid purity. Store Solution II tightly capped when not in use.

* Translation of g (rcf) / rpm see Appendix, page 27

3. Neutralisation of Lysate

Add 350 μ l Solution III to the cleared lysate and gently mix by inverting the tube 6 - 10 times until a flocculent white precipitate is formed. Centrifuge at 10.000 x g* for 10 minutes at room temperature.

4. Load and Bind

Transfer the clear supernatant to a fresh PerfectBind DNA Column in a 2.0 ml Collection Tube. Ensure that the pellet is not disturbed and that no cellular debris are carried over into the column. Centrifuge the column / collection tube assembly for 1 min at 10.000 x g* at room temperature to completely pass lysate through the membrane. Discard the flow-through liquid and keep the Collection Tube for further steps.

5. Wash I (optional)

Place the PerfectBind DNA Column in the 2.0 ml Collection Tube and add 500 μ l PW Plasmid Buffer. Centrifuge for 1 min at 10.000 x g*. Discard flow-through liquid and keep the Collection Tube for further steps.

This wash step with PW Plasmid Buffer leads to an efficient removal of protein contamination and is recommended for sensitive downstream experiments. For standard applications, like cloning screening, you may skip this step.

6. Wash II

Place the PerfectBind DNA Column in the 2.0 ml Collection Tube and add 750 μ l of DNA Wash Buffer completed with ethanol. Centrifuge for 1 min at 10.000 x g*, discard the flow-through and keep the Collection Tube for further steps.

Supplied DNA Wash Buffer must be diluted with absolute ethanol before use. If refrigerated, DNA Wash Buffer must be brought to room temperature before use.

7. Wash III (optional)

Repeat wash step with another 750 μ l DNA Wash Buffer, as described in step 6.

8. Dry (Important! Do not skip this step!)

Place the PerfectBind DNA Column containing your plasmids in the 2.0 ml Collection Tube and centrifuge for 2 min at 10.000 x g* to dry the column matrix. This step is essential to remove ethanol from the column.

* Translation of g (rcf) / rpm see Appendix, page 27

9. Elution

Place the PerfectBind DNA Column into a fresh 1.5 ml microcentrifuge tube. Add 50 - 100 μ l (depending on the desired concentration of final product) Elution Buffer or sterile deionized water directly to the binding matrix in the PerfectBind DNA Column and centrifuge for 1 min at 5.000 x g* to elute DNA.

The first elution represents approximately 75 – 80 % of the bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

Discard the PerfectBind DNA Column and store the eluted plasmid DNA at +4 °C or at -20 °C.

* Translation of g (rcf) / rpm see Appendix, page 27

B. Low Copy-Number Plasmids

1. Bacterial Culture

Inoculate 10 ml growth medium containing the appropriate antibiotic in a 20 - 50 ml culture tube with *E.coli* carrying the desired plasmid and grow at 37 °C with shaking for 12 to 16 h.

Pellet 5 or 10 ml overnight culture by centrifugation for 1 min at 10.000 x g* in several 1.5 ml centrifuge tubes or for 10 min at 5.000 x g* in one 15 ml centrifuge tube. Pour off supernatant and discard.

For Low-Copy Number Plasmids, do not use more than 10 ml overnight culture, otherwise the column may become overloaded. If bigger volumes have to be processed we recommend processing of multiple samples from the same culture. Alternatively, work with the peqGOLD Plasmid Miniprep Kit II (product No. 12-6946 / 12-6945), that allows processing of 25 - 100 ml cultures using the mini-column format and generally yields 20 - 70 µg plasmid DNA.

2. Lysis of Bacteria

Add 250 µl (for 5 ml overnight culture) or 500 µl (for 10 ml overnight culture) Solution I/RNase A to the bacterial pellet, close the tube and completely resuspend cells by vortexing.

Complete resuspension of cell pellet is essential for obtaining high yields.

Add 250 µl (for 5 ml overnight culture) or 500 µl (for 10 ml overnight culture) Solution II and gently mix by inverting and rotating tube 6 - 10 times to obtain a cleared lysate. A 2 min incubation at room temperature may be necessary.

Avoid vigorous mixing as this will shear chromosomal DNA and lower plasmid purity. Store Solution II tightly capped when not in use.

3. Neutralisation of Lysate

Add 350 µl (for 5 ml overnight culture) or 700 µl (for 10 ml overnight culture) Solution III to the cleared lysate and gently mix by inverting the tube 6 - 10 times until a flocculent white precipitate forms. Centrifuge at 10.000 x g* for 10 minutes at room temperature.

4. Load and Bind

Transfer the clear supernatant to a fresh PerfectBind DNA Column in a 2.0 ml Collection Tube. Ensure that the pellet is not disturbed and that no cellular debris are carried over into the column. Centrifuge the column / Collection Tube assembly for 1 min at 10.000 x g* at room temperature to completely pass lysate through the membrane. Discard the flow-through liquid and keep the Collection Tube for further steps.

** Translation of g (rcf) / rpm see Appendix, page 27*

5. Wash I (optional)

Place the PerfectBind DNA Column in the 2.0 ml Collection Tube and add 500 μ l PW Plasmid Buffer. Centrifuge for 1 min at 10.000 x g*. Discard the flow-through liquid and keep the Collection Tube for further steps.

This wash step with PW Plasmid Buffer leads to an efficient removal of protein contamination and is recommended for sensitive downstream experiments. For standard applications, like cloning screening, you may skip this step.

6. Wash II

Place the PerfectBind DNA Column in the 2.0 ml Collection Tube and add 750 μ l of DNA Wash Buffer diluted with ethanol. Centrifuge for 1 min at 10.000 x g*, discard the flow-through and keep the Collection Tube for further steps.

Supplied DNA Wash Buffer must be diluted with absolute ethanol before use. If refrigerated, DNA Wash Buffer must be brought to room temperature before use.

7. Wash III (optional)

Repeat wash step with another 750 μ l DNA Wash Buffer, as described in step 6.

8. Dry (Important! Do not skip this step!)

Place the PerfectBind DNA Column containing your plasmids in the 2.0 ml Collection Tube and centrifuge for 2 min at 10.000 x g* to dry the column matrix. This step is essential to remove ethanol from the column.

* Translation of g (rcf) / rpm see Appendix, page 27

9. Elution

Place the PerfectBind DNA Column into a fresh 1.5 ml microcentrifuge tube. Add 50 - 100 μ l (depending on the desired concentration of final product) Elution Buffer or sterile deionized water directly to the binding matrix in the PerfectBind DNA Column and centrifuge for 1 min at 5.000 x g* to elute DNA.

The first elution represents approximately 75 – 80 % of the bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

Discard the PerfectBind DNA Column and store the eluted plasmid DNA at +4 °C or at –20 °C.

YIELD AND QUALITY OF DNA

Determine the absorption of an appropriate dilution (20- to 50-fold) of the sample at 260 nm and 280 nm. One A_{260} unit corresponds to 50 μ g DNA/ml. The DNA concentration is calculated as follows:

$$\text{DNA conc. } (\mu\text{g/ml}) = \text{Absorption}_{260} \times 50 \times \text{Dilution Factor}$$

High copy number plasmids generally yield up to 25 μ g of DNA from 5 ml culture. The ratio of $A_{260/280}$ is an indicator for nucleic acid purity. A value higher than 1.8 indicates > 90 % purity. Typically, the majority of the eluted DNA is in monomeric supercoiled form.

Alternatively, quantity (as well as quality) can be determined by agarose gel/ethidium bromide electrophoresis by comparison to DNA samples of known concentration.

STORAGE OF DNA

DNA purified by peqGOLD Mini Prep Kits can be stored in Elution Buffer, TE Buffer or sterile, deionized water at –20 °C for years.

* Translation of g (rcf) / rpm see Appendix, page 27

ORDERING INFORMATION

For the isolation of plasmid DNA from bacterial culture:

peqGOLD Plasmid Miniprep Kit I (Classic-Line) <i>(up to 25 µg plasmid DNA)</i>	12-6942-00	5 Preparations
	12-6942-01	50 Preparations
	12-6942-02	200 Preparations
peqGOLD Plasmid Miniprep Kit I (Safety-Line) <i>(up to 25 µg plasmid DNA)</i>	12-6943-00	5 Preparations
	12-6943-01	50 Preparations
	12-6943-02	200 Preparations
peqGOLD Plasmid Miniprep Kit II (Classic-Line) <i>(up to 75 µg plasmid DNA)</i>	12-6945-00	5 Preparations
	12-6945-01	50 Preparations
	12-6945-02	200 Preparations
peqGOLD Plasmid Miniprep Kit II (Safety-Line) <i>(up to 75 µg plasmid DNA)</i>	12-6946-00	5 Preparations
	12-6946-01	50 Preparations
	12-6946-02	200 Preparations
peqGOLD XChange Plasmid Midi Kit <i>(up to 100 µg plasmid DNA)</i>	12-7401-01	10 Preparations
	12-7401-02	20 Preparations
	12-7401-03	40 Preparations
peqGOLD XChange Plasmid Maxi Kit <i>(up to 500 µg plasmid DNA)</i>	12-7402-01	10 Preparations
	12-7402-02	20 Preparations
	12-7402-03	40 Preparations
peqGOLD XChange Plasmid Maxi EF Kit <i>(up to 500 µg plasmid DNA)</i>	12-7404-01	10 Preparations
	12-7404-02	20 Preparations
peqGOLD XChange Purification Kit <i>(for desalination)</i>	12-7403-01	20 Preparations
	12-7403-02	40 Preparations

TROUBLESHOOTING TIPS

Problem	Likely cause	Suggestion
Low DNA yields	Poor cell lysis	<p>Only use LB or YT medium containing appropriate antibiotics</p> <p>Do not use more than 5 ml (HC plasmids) or 10 ml (LC plasmids) culture with the basic protocol</p> <p>Cells may not be suspended good enough prior to addition of Solution II. Vortex thoroughly to completely resuspend the cells</p> <p>Increase incubation time with Solution II to obtain a clear lysate</p> <p>Solution II, if not tightly closed, may need to be replaced. Prepare as follows: 0.1 N NaOH, 1 % SDS</p>
	Bacterial culture overgrown or not fresh	Do not incubate cultures for more than 16 hours at 37 °C and don't use stored cultures
	Low copy-number plasmid used	Such plasmids may yield as little as 0.5 µg DNA from a 5 ml overnight culture. Increase overnight culture volume to 10 ml
No DNA eluted	DNA Wash Buffer Concentrate not diluted with absolute ethanol	Complete DNA Wash Buffer with 100 % ethanol as instructed above
High molecular weight DNA contamination of product	Too strong mixing of cell lysate upon addition of Solution II	Do not vortex or mix aggressively after adding Solution II. Adequate mixing is obtained by simply inverting and rotating tube
Optical densities do not agree with DNA yield estimated by agarose gel-electrophoresis	Trace contaminants eluted from column increase A_{260}	Make sure to wash column as instructed in steps 5 and 7. Alternatively, rely on agarose gel/ethidium bromide electrophoresis for quantitation
RNA visible on agarose gel	RNase A not added to Solution I	Add RNase A to Solution I before use
Plasmid DNA floats out of well while loading agarose gel	Ethanol not completely removed from column following wash steps	Centrifuge column as instructed in step 8 to dry completely

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VWR International GmbH