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INTRODUCTION

The peqGOLD Cycle-Pure Kit is a convenient system for fast and reliable purification of DNA fragments from PCR or other reaction mixes within 15 min. The method uses PerfectBind technology to recover DNA in a range of 50 bp to 40 kb free of oligonucleotides, nucleotides polymerase and other enzymes in yields of 80 % - 90 %. No organic extractions or alcohol precipitations mean safe and rapid processing of multiple samples in parallel.

DNA purified with the peqGOLD Cycle Pure Kit is directly suitable for ligations, PCR sequencing, restriction digestion or various labelling reactions. In addition the kit can be used to purify DNA from any other enzymatic reaction.

peqGOLD Cycle-Pure Kits are available with Safety-Line or Classic-Line columns (Safety-Line, Ordering No. 12-6492-xx; Classic Line, Ordering No. 12-6493-xx). Safety-Line columns can be closed tightly by lids to avoid cross-contamination more effectively. Classic-Line columns do not have lids for a more comfortable handling.

PRINCIPLE

The PCR reaction volume is determined, an equal amount of CP Buffer added to this sample and then applied to a PerfectBind DNA Column. DNA will reversibly bind to the silica matrix and, following two rapid wash steps, where salts, free nucleotides, oligonucleotides and polymerases are removed, the DNA is eluted with deionized water, elution buffer or low salt buffer and ready for further applications.

BENEFITS

- ! Speed - DNA recovery from enzymatic reactions in less than 15 min
- ! Reliability - optimized buffers guarantee pure DNA
- ! Safety - No organic extractions required
- ! Quality - purified DNA suitable for any application

BINDING CAPACITY

Each PerfectBind DNA Column can bind up to ~30 µg DNA.

KIT COMPONENTS

| peqGOLD Cycle Pure Kit | 5 Purifications | 50 Purifications | 200 Purifications |
|--|-----------------|------------------|-------------------|
| Order No. Classic-Line | 12-6493-00 | 12-6493-01 | 12-6493-02 |
| Order No. Safety-Line | 12-6492-00 | 12-6492-01 | 12-6492-02 |
| Components | | | |
| PerfectBind DNA Columns | 5 | 50 | 200 |
| 2 ml Collection Tubes | 5 | 50 | 200 |
| CP Buffer | 1.5 ml | 20 ml | 80 ml |
| CG Wash Buffer | 5 ml | 20 ml | 3 x 20 ml |
| Elution Buffer (10 mM Tris-HCl, pH 9.0) | 1.5 ml | 15 ml | 60 ml |
| Instruction manual | 1 | 1 | 1 |

STORAGE AND STABILITY

All peqGOLD Cycle-Pure Kit components are stable for at least 12 months from the date of purchase when stored at room temperature (22 – 25 °C). Under cool ambient conditions crystals may form in CP Buffer. Simply warm up to 37 °C to dissolve.

BEFORE STARTING

Briefly examine this booklet and become familiar with each step. Prepare all components and have the necessary materials ready before starting. peqGOLD Kits are designed to be simple, fast and reliable if all steps are followed diligently.

! CG Wash Buffer is concentrated and has to be diluted with absolute ethanol as follows:

| | |
|------------|---|
| 12-6492-00 | Add 20 ml 100 % Ethanol to 5 ml CG Wash Buffer |
| 12-6493-00 | |
| 12-6492-01 | Add 80 ml 100 % Ethanol to 20 ml CG Wash Buffer |
| 12-6493-01 | |
| 12-6492-02 | Add 3 x 80 ml 100 % Ethanol to 3 x 20 ml CG Wash Buffer |
| 12-6493-02 | |

! Store diluted CG Wash Buffer at room temperature.

! All steps must be carried out at room temperature.

DANGEROUS COMPONENTS

peqGOLD Cycle-Pure Kit

| Components | Signal word / symbols | Dangerous components | H and P statements |
|--|---|--|---|
| PerfectBind DNA Columns | - | - | - |
| 2 ml Collection Tubes | - | - | - |
| CP Buffer | <p>DANGER</p>  | <p>guanidinium thiocyanate 25-50%, propan-2-ol 20-<25%, polyethylene glycol octylphenol ether 0.3-<1%, CAS: 593-84-0, 67-63-0, 9002-93-1</p> | <p>H225, H302+H312+H332, H314, H336, H412, P101, P102, P103, P210, P241, P303+P361+P353, P305+P351+P338, P405, P501</p> |
| CG Wash Buffer | - | - | - |
| Elution Buffer (10 mM Tris-HCl, pH 9.0) | - | - | - |

| H and P statements | Descriptions |
|---------------------------|--|
| H225 | Highly flammable liquid and vapour. |
| H302+H312+H332 | Harmful if swallowed, in contact with skin or if inhaled. |
| H314 | Causes severe skin burns and eye damage. |
| H336 | May cause drowsiness or dizziness. |
| H412 | Harmful to aquatic life with long lasting effects. |
| P101 | If medical advice is needed, have product container or label at hand. |
| P102 | Keep out of reach of children. |
| P103 | Read label before use. |
| P210 | Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. |
| P241 | Use explosion-proof electrical/ventilating/lighting/.../ equipment. |
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P405 | Store locked up. |
| P501 | Dispose of contents/container to ... |

PEQGOLD CYCLE-PURE PROTOCOL

Materials to be supplied by user:

- ! Sterile 1.5 ml centrifuge tubes
- ! Sterile deionized water (optional)
- ! 100 % ethanol

1. Agarose gel electrophoresis

Perform agarose gel/ethidium bromide electrophoresis to analyze PCR product.

2. Load and Bind

Determine the volume of the PCR reaction, transfer to a clean 1.5 ml microfuge tube, add an equal volume of CP Buffer and vortex thoroughly to mix.

For PCR products < 200 bp add 3 volumes of CP Buffer.

Apply the sample to a PerfectBind DNA Column assembled in a clean 2 ml Collection Tube (provided) and centrifuge in a microcentrifuge at 10.000 x g* for 1 min at room temperature. Discard the liquid and reuse collection tube in the following steps.

3. Wash

Wash the PerfectBind DNA Column by adding 750 µl of CG Wash Buffer diluted with absolute ethanol. Centrifuge at 10.000 x g* for 1 min at room temperature. Discard liquid and repeat this step.

4. Dry (Important, do not skip this step!)

Discard liquid and centrifuge the empty column for 1 min 10.000 x g* to dry the column matrix. This step is essential for good DNA yields.

5. Elution

Place PerfectBind DNA Column into a clean 1.5 ml microcentrifuge tube. Add 30 – 50 µl Elution Buffer (depending on desired concentration of final product) or sterile deionized water directly onto the column matrix and centrifuge 1 min at 5.000 x g* to elute DNA. This represents approximately 80 – 90 % of bound DNA. An optional second elution will improve yield up to 95 %, though at a lower concentration.

For fragments > 500 bp you usually get more than 80 % yield. For fragments of 50 bp – 500 bp you usually get 55 % - 80 %.

* Translation g (rcf) / rpm see Appendix, p. 19

CONCENTRATION OF DNA

If required DNA fragments purified with peqGOLD Cycle-Pure Kits can be further concentrated. Therefore add NaCl up to the final concentration of 0.1M and then add 2 volumes of 100 % ethanol. Vortex thoroughly and centrifuge at 10.000 x g* for 15 min. Discard the supernatant, add 700 µl 70 % Ethanol and centrifuge again at 10.000 x g* for 2 min. Discard the supernatant and air-dry the pellet. Solve the dried DNA in 20 µl sterile deionized water, Elution or other low salt buffer.

YIELD AND QUALITY OF DNA

Determine the absorption of an appropriate dilution (20- to 50-fold) of the sample at 260 nm and 280 nm. One A_{260} -unit corresponds to 50 µg DNA/ml. The DNA concentration is calculated as follows:

$$\text{DNA conc. } (\mu\text{g/ml}) = \text{Absorption}_{260} \times 50 \times \text{Dilution Factor}$$

The ratio of $A_{260/280}$ is an indicator for nucleic acid purity. A value 1.8 – 2.0 indicates a purity of 90 % - 100 %.

Alternatively, quantity (as well as quality) can sometimes best be determined by agarose gel/ethidium bromide electrophoresis by comparison to DNA samples of known concentrations.

DNA purified with the peqGOLD Cycle-Pure Kit can be stored in Elution Buffer or sterile, deionized water at –20 °C for years.

* Translation g (rcf) / rpm see Appendix, p. 19

SHORT PROTOCOL FOR EXPERIENCED USERS

1. Determine volume of PCR reaction.

2. Add equal volume of CP Buffer to PCR reaction.

3. Apply 750 μ l of PCR reaction/CP Buffer solution to the PerfectBind DNA Column assembled in a 2 ml Collection Tube.

4. Wash column twice with 750 μ l CG Wash Buffer.

5. Dry column by centrifugation.

6. Place column in a fresh 1.5 ml tube and elute DNA with 30 – 50 μ l of Elution Buffer.

ORDERING INFORMATION

For extraction of DNA from agarose gels and PCR reactions:

| | | |
|---|------------|------------------|
| peqGOLD Gel Extraction Kit | 12-2501-00 | 5 Preparations |
| (Classic-Line) | 12-2501-01 | 50 Preparations |
| <i>(DNA from agarose gels)</i> | 12-2501-02 | 200 Preparations |
| <hr/> | | |
| peqGOLD Gel Extraction Kit | 12-2500-00 | 5 Preparations |
| (Safety-Line) | 12-2500-01 | 50 Preparations |
| <i>(DNA from agarose gels)</i> | 12-2500-02 | 200 Preparations |
| <hr/> | | |
| peqGOLD Cycle-Pure Kit | 12-6493-00 | 5 Preparations |
| (Classic-Line) | 12-6493-01 | 50 Preparations |
| <i>(DNA from PCR reactions)</i> | 12-6493-02 | 200 Preparations |
| <hr/> | | |
| peqGOLD Cycle-Pure Kit | 12-6492-00 | 5 Preparations |
| (Safety-Line) | 12-6492-01 | 50 Preparations |
| <i>(DNA from PCR reactions)</i> | 12-6492-02 | 200 Preparations |
| <hr/> | | |
| peqGOLD MicroSpin Gel Extraction Kit | 12-6294-00 | 5 Preparations |
| (Safety-Line) | 12-6294-01 | 50 Preparations |
| <i>(DNA from agarose gels)</i> | 12-6294-02 | 200 Preparations |
| <hr/> | | |
| peqGOLD MicroSpin Cycle-Pure Kit | 12-6293-00 | 5 Preparations |
| (Safety-Line) | 12-6293-01 | 50 Preparations |
| <i>(DNA from PCR reactions)</i> | 12-6293-02 | 200 Preparations |

TROUBLESHOOTING TIPS

| Problem | Likely cause | Suggestion |
|--|---|---|
| Low DNA yields | Too little CP Buffer added to sample | Add CP Buffer as indicated in the manual |
| No DNA eluted | CG Wash Buffer Concentrate not diluted with absolute ethanol | Complete CG Wash Buffer Concentrate as instructed in the manual with 100 % Ethanol |
| | Too little CP Buffer added to sample | Add CP Buffer as indicated in the manual |
| Optical densities do not agree with DNA yield on agarose gel | Trace contaminants eluted from column increase A_{260} | Wash column twice as instructed in step 3. Alternatively, perform agarose gel/ethidium bromide electrophoresis for quantification |
| DNA floats out of well while loading agarose gel | Ethanol not completely removed from column following wash steps | Centrifuge column as instructed in step 4 to dry before proceeding to elution step |

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VWR International GmbH