BrdU Cell Proliferation Assay Kit #6813

Protocol

BrdU Cell Proliferation Assay Kit Protocol

A. Reagent Preparation

- 1. Prepare 1X Wash Buffer by diluting 20X Wash Buffer (included in each BrdU ELISA Kit) in purified water.
- 2. Prepare 1X detection antibody solution by diluting BrdU Detection Antibody 1:100 with Detection Antibody Diluent (green).
- 3. Prepare 1X HRP-conjugated secondary antibody solution by diluting Anti-mouse IgG, HRP-linked Antibody 1:100 with HRP-linked Antibody Diluent (red).
- 4. Prepare 10X BrdU solution by diluting BrdU 1:100 with cell culture medium.

B. BrdU Incorporation

- 1. Plate cells in 96-well plate and incubate with respective test substance. Typical seed cell number is 2500– 100000 cells/well depending on cell growth rate. Typical incubation time is 1–72 hr.
- 2. Add prepared 10X BrdU solution to plate wells, for a final 1X concentration. (Example: For 100 μ l medium in the plate, add 10 μ l of 10X BrdU solution per well.)
- 3. Place cells in incubator. Typical incubation time is 1–24 hr.
- 4. Remove medium. For suspension cells, centrifuge the plate at 300 g for 10 min, then remove medium.

C. BrdU Assay

- 1. Add 100 µl/well of the Fixing/Denaturing Solution, keep the plate at room temperature for 30 min. Remove solution.
- 2. Add 100 µl/well prepared 1X detection antibody solution, keep plate at room temperature for 1 hr. Remove solution and wash plate 3 times with 1X Wash Buffer.
- 3. Add 100 µl/well prepared 1X HRP-conjugated secondary antibody solution, keep plate at room temperature for 30 min. Remove the solution and wash plate 3 times with 1X Wash Buffer.
- 4. Add 100 µl TMB Substrate (#7004 (/product/productDetail.jsp?productId=7004)).
- 5. Incubate for 30 min at room temperature.

NOTE: Watch the color change as it may be necessary to stop the reaction prior to the standard development time of 30 min.

- 6. Add 100 µl STOP Solution (#7002 (/product/productDetail.jsp?productId=7002)).
- 7. Read absorbance at 450 nm (For optimal readings, read the plate within 30 min of adding STOP Solution).

posted October 2012

Λ