

# RIPA Buffer

[Learn More](#)

89900 89901

1782.1

Number	Description
89900 <a href="#">Buy Now</a>	RIPA Buffer, 100mL
89901 <a href="#">Buy Now</a>	RIPA Buffer, 250mL

Contents: 25mM Tris•HCl pH 7.6, 150mM NaCl, 1% NP-40, 1% sodium deoxycholate, 0.1% SDS

**Storage:** Upon receipt store at 4°C. Product shipped at ambient temperature.

## Introduction

The Thermo Scientific RIPA buffer is one of the most reliable buffers used to lyse cultured mammalian cells from both plated cells and cells pelleted from suspension cultures. This buffer enables protein extraction from cytoplasmic, membrane and nuclear proteins and is compatible with many applications, including reporter assays, protein assays, immunoassays and protein purification.

## Important Product Information

- RIPA Buffer does not contain protease or phosphatase inhibitors. If desired, add protease inhibitors, such as Thermo Scientific Halt Protease Inhibitor Cocktail (Product No. 78410) and Halt™ Phosphatase Inhibitor Cocktail (Product No. 78420) to the reagent to prevent proteolysis and maintain phosphorylation status of proteins. Add protease and phosphatase inhibitors immediately before use.
- Use 1mL of cold RIPA Buffer for every  $5 \times 10^6$  of HeLa or A431 cells (~20μL of packed cells, which is equivalent to ~40mg of cells). To obtain concentrated protein extracts, directly lyse cells on plate and use less buffer.
- Some protein kinases and other enzymes may be sensitive to the components of the RIPA Buffer, resulting in their decreased activity. In such cases, prepare a RIPA buffer that does not contain sodium deoxycholate and SDS.
- RIPA Buffer is compatible with the Thermo Scientific Pierce BCA Protein Assay Kit (Product No 23225).

## Procedure for Lysis of Monolayer-cultured Mammalian Cells

**Note:** If desired, add protease and phosphatase inhibitors to the RIPA Buffer immediately before use.

1. Carefully remove (decant) culture medium from adherent cells.
2. Wash cells twice with cold PBS.
3. Add cold RIPA Buffer to the cells. Use 1mL of buffer per 75cm<sup>2</sup> flask containing  $5 \times 10^6$  HeLa or A431 cells. Keep on ice for 5 minutes, swirling the plate occasionally for uniform spreading.
4. Gather the lysate to one side using a cell scraper, collect the lysate and transfer to a microcentrifuge tube. Centrifuge samples at ~14,000 × g for 15 minutes to pellet the cell debris.

**Note:** To increase yields, sonicate the pellet for 30 seconds with 50% pulse.

5. Transfer supernatant to a new tube for further analysis.

## Procedure for Lysis of Suspension-cultured Mammalian Cells

**Note:** If desired, add protease and phosphatase inhibitors to the RIPA Buffer immediately before use.

1. Pellet the cells by centrifugation at  $2500 \times g$  for 5 minutes. Discard the supernatant.
2. Wash cells twice in cold PBS. Pellet cells by centrifugation at  $2500 \times g$  for 5 minutes.
3. Add RIPA Buffer to the cell pellet. Use 1mL of RIPA buffer for 40mg ( $\sim 5 \times 10^6$  of HeLa cells) of wet cell pellet. Pipette the mixture up and down to suspend the pellet.

**Note:** To increase yields, sonicate the pellet for 30 seconds with 50% pulse.

4. Shake mixture gently for 15 minutes on ice. Centrifuge mixture at  $\sim 14,000 \times g$  for 15 minutes to pellet the cell debris.
5. Transfer supernatant to a new tube for further analysis.

## Troubleshooting

Problem	Possible Cause	Solution
Low total protein yield	Some cells are more resistant to lysis than others	Make sure the cell pellet is thoroughly suspended in RIPA Buffer and incubate for longer with occasional swirling – sonicate the pellet to increase yield
Low concentration of proteins	Excess buffer used	Use less buffer (e.g., 0.25-0.5mL per 75cm <sup>2</sup> flask containing $5 \times 10^6$ cells) – use a sufficient amount to cover the entire plate
Proteolysis	No protease inhibitors added	Add Halt Protease Inhibitor Cocktail to the buffer before use
Low phosphorylation of proteins	Phosphatase activity	Add Halt Phosphatase Inhibitor Cocktail to the buffer before use
	Protein is non-phosphorylated or poorly phosphorylated	None

## Related Thermo Scientific Products

78410	<b>Halt Protease Inhibitor Cocktail Kit</b>
78420	<b>Halt Phosphatase Inhibitor Cocktail, 1mL</b>
78248	<b>B-PER<sup>®</sup> Bacterial Protein Extraction Reagent, 500mL</b>
78990	<b>Y-PER<sup>®</sup> Yeast Protein Extraction Reagent, 500mL</b>
89826	<b>Mem-PER<sup>®</sup> Membrane Protein Extraction Reagent Kit</b>
78833	<b>NE-PER<sup>®</sup> Nuclear and Cytoplasmic Extraction Kit</b>
23227	<b>Pierce<sup>®</sup> BCA Protein Assay Kit</b>
26148	<b>Pierce Direct IP Kit</b>
34080	<b>SuperSignal<sup>®</sup> West Pico Chemiluminescent Substrate, 500mL</b>
34076	<b>SuperSignal<sup>®</sup> West Dura Extended Duration Substrate, 200mL</b>

## General References

- Cao, F., *et al.* (2005). Identification of an essential molecular contact point on the duck hepatitis B virus reverse transcriptase. *J Virol* **79**(16):10164-70.
- Pfeffer, K.I. and Flugel R.M. (2005). Molecular characterization of proteolytic processing of the gap proteins of human spumaretrovirus. *Methods in Mol Biol* **304**:435-44.
- Sefton, B.M. (2005). Labeling cultured cells with 32Pi and preparing cell lysates for immunoprecipitation. Unit 18.2. F. M. Ausubel, R. Brent, R.E. Kingston, D.D. Moore, J.G. Seidman, J.A. Smith, and K. Struhl (eds.) *Current Protocols in Molecular Biology*. John Wiley & Sons, Inc.

---

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.