

# Microbial Growth

## Related Links

[VIEW PRODUCT LIST](#)

## Categories

- [Antibiotics](#)

## Microbial Growth Protocols

Find more protocols and selection guides in the [Molecular Biology Guide](#).

### Culture of *E. coli*

#### Materials required:

Reagents	Equipment
Pure <i>E. coli</i> culture	Cabinet incubator (37°C)
LB – Miller	Shaker incubator (37°C)
LB – Miller Plates	Culture dishes, sterile
	Culture flasks, sterile

#### Suspension culture of *E. coli*

Prepare LB medium by weighing appropriate powder medium and adding to water in a sterile flask
↓
Autoclave the broth and cool to room temperature (Alternatively ready-to-use LB medium may be used)
↓
In a laminar flow chamber, transfer approximately 1 mL of overnight <i>E. coli</i> culture to the flask
↓
Seal the mouth of the flask with sterile cotton (non-absorbant) plugs; ensure the flask is not tightly sealed
↓
Incubate overnight at 37°C with continuous shaking

#### Plating *E. coli*

Prepare LB agar by weighing appropriate powder medium, agar and water to a sterile flask		
↓		
Autoclave the medium and cool just enough to be able to handle the flask		
↓		
In a laminar flow chamber, pour approximately 25-30 mL of the LB agar into sterile plates		
↓		
Allow the plates to set in the laminar flow chamber with lids slightly opened (The plates may also be stored inverted at 4°C for future use) (Alternatively, ready-to-use LB agar plates may be used)		
↓	↓	↓
Lightly scratch the surface of frozen <i>E. coli</i> glycerol stock with a sterile inoculating loop	Pick up <i>E. coli</i> colony from a plate with culture with a sterile inoculating loop	Add 10-100 µL of <i>E. coli</i> suspension culture and add one the LB agar plate
↓	↓	↓
Streak the loop across the LB agar plate		Spread the culture all over the plate using a sterile glass spreader

↓  
Invert and incubate the plates overnight at 37°C  
↓

## Advanced Protocols

### Plating Bacteriophage M13

#### Materials required:



Reagents	Equipment
Pure <i>E. coli</i> culture	Cabinet incubator (37°C)
<b>LB agar EZMix™ (Product No. L7533)</b>	Water Bath (47°C)
Appropriate selection antibiotic	Culture dishes
<b>IPTG (Isopropyl-β-D-thiogalactoside) (Product No. I6758)</b>	15 ml polypropylene culture tubes (sterile)
<b>X-gal (Product No. B9146)</b>	
Magnesium chloride	






Prepare a pure culture of <i>E. coli</i> by inoculating a single colony into 5 mL of LB or YT medium
↓
Prepare ten-fold serial dilutions of M13 bacteriophage stock in LB or YT medium
↓
Prepare LB agar or YT medium supplemented with 5mM MgCl <sub>2</sub> in sterile test tubes; equilibrate at 47°C; add appropriate amounts of X-gal and IPTG solutions
↓
Add serially diluted M13 bacteriophage stocks 100 μL of pure bacterial culture; mix by gentle vortexing
↓
Pour the LB agar or YT medium with X-gal and IPTG to tubes containing infected bacteria; mix by gentle vortexing
↓
Transfer the contents to plate and swirl for even distribution of infected bacteria
↓
Allow the plates to set; invert and incubate the plates at 37°C
↓
Pale blue plaques of M13 bacteriophage appear on a lawn of bacterial growth

The following are the bacterial broth components offered by Sigma-Aldrich:

Bacterial Broth Components	Name	Format
<b>Y0875</b>	<b>Select Yeast Extract</b>	Powder
<b>N4767</b>	<b>N-Z-Amine® EKC</b>	
<b>Y1626</b>	<b>Select Yeast Extract, EZMix™ powder</b>	EZMix™
<b>T2559</b>	<b>EZMix™ Tryptone</b>	

## Materials


Product #	Image	Description	Molecular Formula	Add to Cart
<b>Y2627</b>		2xYT Medium EZMix™ Powder microbial growth medium		<a href="#">pricing</a>
<b>B9146</b>		5-Bromo-4-chloro-3-indolyl β-D-galactopyranoside ≥98%	C <sub>14</sub> H <sub>15</sub> BrClNO <sub>6</sub>	<a href="#">pricing</a>
<b>C4464</b>		EZMix™ N-Z-Amine® A Casein enzymatic hydrolysate, from bovine milk		<a href="#">pricing</a>
<b>T2559</b>		EZMix™ Tryptone Pancreatic digest of casein		<a href="#">pricing</a>
<b>H8032</b>		Hanahan's Broth (SOB Medium) powder microbial growth medium		<a href="#">pricing</a>
<b>I6758</b>		IPTG ≥99% (TLC), ≤0.1% Dioxane	C <sub>9</sub> H <sub>18</sub> O <sub>5</sub> S	<a href="#">pricing</a>
<b>L3022</b>		LB Broth (Lennox) Highly-referenced microbial growth powder medium, low salt, suitable for salt-sensitive <i>E. coli</i> culture.		<a href="#">pricing</a>
<b>L7275</b>		LB Broth (Lennox) Tablet microbial growth medium		<a href="#">pricing</a>

Product #	Image	Description	Molecular Formula	Add to Cart
L7658		LB Broth (Lennox) EZMix™ powder microbial growth medium		<a href="#">pricing</a>
L3397		LB Broth (Luria low salt) Powder microbial growth medium		<a href="#">pricing</a>
L3522		LB Broth (Miller) Highly-referenced nutrient-rich microbial growth powder medium, suitable for regular <i>E.coli</i> culture		<a href="#">pricing</a>
L2542		LB Broth (Miller) Liquid microbial growth medium		<a href="#">pricing</a>
L2897		LB Broth with agar (Lennox) Highly-referenced microbial growth powder medium with Agar, low salt, suitable for salt-sensitive <i>E. coli</i> culture.		<a href="#">pricing</a>
L7025		LB Broth with agar (Lennox) Tablet microbial growth medium		<a href="#">pricing</a>
L7533		LB Broth with agar (Lennox) EZMix™ powder microbial growth medium		<a href="#">pricing</a>
L3272		LB Broth with agar (Luria low salt) Powder microbial growth medium		<a href="#">pricing</a>
L3147		LB Broth with agar (Miller) Highly-referenced nutrient-rich microbial growth powder medium with Agar, suitable for regular <i>E.coli</i> culture.		<a href="#">pricing</a>
M6030		M9, Minimal Salts, 5X powder, minimal microbial growth medium		<a href="#">pricing</a>
N3518		NZ Amine® Broth Lambda bacteriophage growth medium		<a href="#">pricing</a>
N4767		N-Z-Amine® EKC Casein enzymatic hydrolysate, from bovine milk		<a href="#">pricing</a>
N3643		NZCYM Broth Powder microbial growth medium		<a href="#">pricing</a>
N6905		NZCYM Broth EZMix™ powder microbial growth medium		<a href="#">pricing</a>
Y0875		Select Yeast Extract for use in microbial growth medium		<a href="#">pricing</a>
Y1626		Select Yeast Extract, EZMix™ powder for use in microbial growth medium		<a href="#">pricing</a>
S1797		SOC Medium For use in transformation		<a href="#">pricing</a>
T9179		Terrific Broth EZMix™ powder microbial growth medium		<a href="#">pricing</a>
T5574		Terrific Broth Liquid microbial growth medium		<a href="#">pricing</a>
Y2377		2x YT medium Powder microbial growth medium		<a href="#">pricing</a>


## References

- Sambrook, J., Fritsch, E.F., and Maniatis, T., Molecular cloning: a laboratory manual. New York: Cold spring harbor laboratory press, 1989.
- Parija, SC., Textbook of Microbiology and Immunology. Elsevier India, 2009.
- <http://www.scienceprofonline.org/microbiology/types-culture-media-for-growing-bacteria.html>

**Web Toolbox**  
Calculators | Selectors | Explorers

**General Help**  
 [Email Customer Support](#)

**Structure Search**  
JME Editor  
Courtesy of Peter Ertl

**Ask a scientist**  
 [Email Technical Service](#)

### Service & Support

CUSTOMER SUPPORT  
TECHNICAL SERVICE  
WEB HELP DESK  
SDS  
C OF A

### Ordering

CUSTOM PRODUCTS  
ECOMMERCE SOLUTIONS  
ORDER CENTER  
PRODUCTS  
TERMS & CONDITIONS OF SALE

### Corporate

BUSINESS DEVELOPMENT  
WORLDWIDE OFFICES  
ABOUT US  
SITE MAP  
CAREERS  
EVENTS  
PROGRAMS  
REACH REGULATIONS  
CONTACT US  
EMAIL SUBSCRIPTION CENTER  
TOOL BOX

