TaqMan® Pathogen Detection Kit

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Sep, 2006

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The screenshot shows a form titled 'Enter Plate Information'. It contains several input fields: 'Plate Name' (with 'Plate_02142006_1' entered), 'Barcode', 'Setup Operator', and 'Comments' (a large text area). A note next to the Barcode field says 'To scan the barcode, click the Barcode field, then press the trigger on the scanner.'

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Select Target(s)
Select the target pathogens to test for, then enter the number of samples, replicates, and positive and negative controls for each required for each target pathogen. Click Next to continue.

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Enter Sample Names
Enter or import sample names and information for each target. To change the number of samples or replicates for a target, click

Plate Name: Plate_02142006_1

Target/Sample Information

To add additional targets or samples, click Back to return to the previous page.

Target	Sample Name	Replicates	User Defined #1	User Defined #2	User Defined #3
2	E. coli O157:H7	No.2			
3	E. coli O157:H7	No.3			
4	E. coli O157:H7	No.4			
5	E. coli O157:H7	No.5			
6	S. enterica	No.1			
7	S. enterica	No.2			
8	S. enterica	No.3			
9	S. enterica	No.4			
10	S. enterica	No.5			


Fill Down Delete Selected

Target/Control Information


Target	Control	Replicates	
1	E. coli O157:H7	Negative	3
2	S. enterica	Negative	3

This illustration shows the number of wells containing samples
Wells Used: 36 / 96

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


Review Reagent and Sample Information
Review the information listed for reagents and samples. Ensure that you have the indicated volume of all items, check the box at the bottom

Plate Information
Plate Name: Plate_02142006_1 Barcode: Created: 2006/02/14 15:58

Reagents

Pipette Overage: 5 %

	Name	Kit Part Number	Volume
1	E. coli O157:H7	4366100	
2	EMM	kit component	283.50 uL
3	TAM	kit component	56.70 uL
4	Negative Control	kit component	36.00 uL
5	S. enterica	4366104	
6	EMM	kit component	283.50 uL
7	TAM	kit component	56.70 uL
8	Negative Control	kit component	36.00 uL

Samples

	Name	Volume
1	No.1	72 uL
2	No.2	72 uL
3	No.3	72 uL
4	No.4	72 uL
5	No.5	72 uL

10

Ensure all reagents and samples are available

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
Prepare Premix Solutions
Follow the directions below to prepare the needed premix solutions, check the box at the bottom of the page, then click Next.

Plate Information
Plate Name: Plate_02142006_1 Barcode: Created: 2006/02/14 15:58:11

Premix Directions

1. Thaw all reagents completely, mix (by gently aspirating/dispensing a few times or by vortexing briefly), then centrifuge briefly.
2. Label a tube with target name.
3. Combine the volumes listed of 2x Environmental Master Mix (EMM) and 10x Target Assay Mix (TAM). Use a new tip when pipetting EMM and TAM.
4. Mix reagents (by aspirating/dispensing a few times or by vortexing briefly).
5. Close the tube.
6. Repeat steps 1 through 5 for each target assay.
7. When all premix solutions are prepared, check the box at the bottom of the page, then click Next.

	Target Name	2x EMM Vol	10x TAM Vol
1	E. coli O157:H7	283.50 uL	56.70 uL
2	S. enterica	283.50 uL	56.70 uL

Select to confirm that this step has been completed

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Plate Information
Plate Name: Plate_02142006_1 Barcode: Created: 2006/02/14 15:58:11

Pipette Premix
1. Click a Target Name in the table at the right to highlight the wells reserved for the premix solution for the selected target.
2. Pipette the indicated volume of premix solution for the selected target into the highlighted wells.
3. (optional) When all wells for the selected target are pipetted check the Pipetted box in the table at the right.
4. Change the pipette tip, then repeat steps 1 through 3 for the next target.
5. When all premix solutions are pipetted, check the box at the bottom of the page, then click Next.

Pipetted	Target Name	Volume/W
<input checked="" type="checkbox"/>	E. coli O157:H7	18.00 μ L
<input checked="" type="checkbox"/>	S. enterica	18.00 μ L

Select to confirm that this step has been completed

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Pipette Unknown Samples
Follow the directions below to pipette unknown samples, check the box at the bottom of the page, then click Next.

Plate Information
Plate Name: Plate_02142006_1 Barcode: Created:

Pipette Unknowns
Target (1 of 2): E. coli O157:H7
The unknown sample wells reserved for the target above are highlighted below.

1. Click a Well Number in the table at right, then pipette the indicated volume of an unknown into the indicated well.
2. Mix by gently aspirating and dispensing a few times. IMPORTANT! Mix very gently with the pipette tip at the bottom of the tube to minimize aerosol formation and cross-contamination.
3. (optional) Check the Pipetted box.
4. Change the pipette tip, then repeat steps 1 through 4 for the next unknown.
5. When all unknown samples for the target are pipetted, check the box at the bottom of the page, then click Next.

Pipetted	Well Number	Sample No
<input checked="" type="checkbox"/>	A1	No.1
<input checked="" type="checkbox"/>	A2	No.1
<input checked="" type="checkbox"/>	A3	No.1
<input checked="" type="checkbox"/>	A4	No.2
<input checked="" type="checkbox"/>	A5	No.2
<input checked="" type="checkbox"/>	A6	No.2
<input type="checkbox"/>	A7	No.3
<input type="checkbox"/>	A8	No.3
<input type="checkbox"/>	A9	No.3

Select to confirm that this step has been completed

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Close Pipetted Tubes Before Pipetting Positive Controls
Follow the directions below to close tubes, check the box at the bottom of the page, then click Next.

Plate Information
Plate Name: Plate_07052006_1 Barcode: Created: 2006/07/05 11:31

Close Tubes
Close all unknown sample and negative control tubes.
IMPORTANT! If you are using an optical cover on a plate, skip this step, and be extremely careful to avoid contamination when pipetting positive controls.

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Pipette Positive Controls
Follow the directions below to pipette positive controls, check the box at the bottom of the page, then click Next.

Plate Information
Plate Name: Plate_07052006_1 Barcode:

Pipette Positive Controls
IMPORTANT! When pipetting positive control, be extremely careful. Contamination of unknown wells with positive control will affect results.

1. Click a Well Number in the table at the right, then pipette the indicated volume of the positive control for the target into the indicated well.
2. Mix by gently aspirating and dispensing a few times. **IMPORTANT!** Mix very gently with the pipette tip at the bottom of the tube to minimize aerosol formation and cross-contamination.
3. (optional) Check the Pipetted box.
4. Change the pipette tip, then repeat steps 1 through 3 for the next positive control well.
5. When all positive controls are pipetted, check the check box at the bottom of the page, then click Next.

Pipetted	Well Number
<input checked="" type="checkbox"/>	B11
<input type="checkbox"/>	B12
<input type="checkbox"/>	E11
<input type="checkbox"/>	E12

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Close Tubes
Follow the directions below to close tubes, check the box at the bottom of the page, then click Next.

Plate Information
Plate Name: Plate_07052006_1 Barcode: Created: 2006/07

Close Tubes

1. Close all tubes or apply an optical cover to the plate.
2. Make sure reagents are in the bottom of the wells. If available, use a centrifuge with a plate adapter to briefly centrifuge the plate (<1500g).
3. Check the check box at the bottom of the page, then click Next.

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Confirm Plate Map
Review the plate map, then click Next to confirm.

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Select the Next Step
Select your next step.

Plate Information
Plate Name: Plate_02142006_1 Barcode:

Next Step Options

- Run plate**
Select this option to run a plate.
- Print pipetting procedures**
Select this option to print the pipetting procedure you just performed.
- Pipette another plate**
Select this option to pipette another plate.

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Select a Plate to Run
Select a plate, review the plate details, and add information as needed.
Load the plate in the instrument (click Help for information), then click Next.

Select a Plate
Date Range: All Dates
(mm/dd/yyyy) (mm/dd/yyyy)
From: 01/01/1980 To: 02/14/2006

Plate Name /#	Last Modified
Plate_01272006_1	01/31/2006
Plate_01312006_1	01/31/2006
Plate_01312006_2	01/31/2006
Plate_02102006_1	02/10/2006
Plate_02102006_2	02/10/2006
Plate_02102006_3	02/10/2006
Plate_02142006_1	02/14/2006

Plate Details
Enter Run Operator name (optional) and Comment (optional).

Plate Name: Plate_02142006_1
Barcode:
Setup Operator:
Plate Setup: 2006/02/14 15:58:11
Pipette Operator:
Run Operator:
Instrument Type: 7500 System
Software Version: 1.0
Targets on Plate: E. coli O157:H7; S. enterica
Comments:

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 **Run a Plate**
Change the Number of PCR Cycles as needed, then click Start.

Plate Information
Plate Name: Plate_02142006_1 Barcode:

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
View Results
Select a plate, then click the tabs below to view results for wells or samples. To print or export results, click the Report tab.

Plate: Complete Runs

Range: All Dates

(mm/dd/yyyy) (mm/dd/yyyy)
01/01/1980 To: 02/14/2006

Name / #	Crt Date	Run Date
ple_Ecoli_Sal	02/14/20...	07/29/2005
02142006_1	02/14/20...	02/14/2006
	01/27/20...	07/29/2005

Results by Well | Results By Sample | Plate Information | Report

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The screenshot displays the '7900 pathogen picture.sds - Absolute Quantification' software interface. The main area is a 4x10 grid of wells (A-D, 1-10). Each well contains a colored box with a 'U' and a detector name. For example, well A1 has 'U IPC' and 'U Liste'. The right panel shows the configuration for the selected wells. The 'Sample Name' field is set to '* Mixed *'. Below it is a table with the following columns: Use, Detector, Reporter, Task, Quantity, and Color.

Use	Detector	Reporter	Task	Quantity	Color
<input checked="" type="checkbox"/>	IPC	VIC	Unkno...	0	Blue
<input checked="" type="checkbox"/>	Listeria	FAM	Unkno...	0	Red
<input type="checkbox"/>	Salmonella	FAM		0	Yellow

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The PCR process and 5' nuclease process are covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd.

Applied Biosystems, PrepMan and VIC are registered trademarks and AB (Design), Applera, FAM and are trademarks of Applera Corporation or its subsidiaries in the US and/or certain other countries.

Rapid Finder is a trademark of Applied Biosystems/MDS SCIEX, which is a joint venture between Applera Corporation and MDS Inc.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.

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© 2006 Applied Biosystems Japan Ltd. All rights reserved.

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