remel Candida Albicans Test Kit

REF R30851401.....50 Tests/Kit

1. INTENDED USE

Remel Candida Albicans Test Kit is recommended for use in qualitative procedures to aid in the identification of *Candida albicans* by detection of the enzymes L-proline aminopeptidase and β -galactos-aminidase in yeast cells cultured on laboratory media.

2. SUMMARY AND EXPLANATION

C. albicans is the species of yeast most frequently isolated from clinical specimens. Historically, *C. albicans* has been presumptively identified using the germ tube test.¹ Studies, however, have demonstrated that up to 5% of *C. albicans* clinical isolates are germ tube negative. Additionally, early pseudohyphae of *Candida tropicalis* may be confused with germ tube formation.² In 1987, Perry and Miller reported the use of an umbelliferyl-labeled compound (N-acetyl- β -D-galactosaminide) to differentiate *C. albicans* and *C. tropicalis*. This test was later modified to include the detection of L-proline aminopeptidase to increase the test specificity. Also, the substitution of a nitrophenyl-labeled compound provided a colorimetric rather than fluorometric endpoint for the detection of β -galactosaminidase.³ While other species of yeasts were found to be positive for L-proline aminopeptidase or β -galactosaminidase, *C. albicans* was positive for both enzymes. This substrate combination provides a rapid, non-microscopic method to aid in the identification of *C. albicans*.⁴

3. PRINCIPLE

The substrates, p-nitrophenyl-N-acetyl- β -D-galactosaminide (NGL) and L-proline- β -naphthylamide (PRO) are impregnated on each Candida Albicans Test Kit disk. The hydrolysis of NGL by the enzyme β -galactosaminidase is indicated by the development of a yellow color due to the liberation of p-nitrophenol. The addition of Sodium Hydroxide, 0.3% enhances the reaction by creating alkaline conditions. If the PRO substrate has been hydrolyzed by L-proline aminopeptidase, a red color develops due to the reaction of β -naphthylamine with p-dimethylaminocinnamaldehyde.

4. REAGENTS AND MATERIALS SUPPLIED

- Candida Albicans Test Kit Disks (1 vial of 50)
Reactive Ingredient: p-nitrophenol-N-acetyl- β -D-galactosaminide (NGL) and L-proline- β -naphthylamide (PRO)
- Color Developer (1 bottle, 4.5 ml)
Reactive Ingredient: 0.01% p-dimethylaminocinnamaldehyde
- Sodium Hydroxide, 0.3% (1 bottle, 4.5 ml)
- Plastic Tubes with Caps (50)
- Applicator Sticks (1 vial of 50)
- Instructions for Use (IFU)

5. CONTENTS SYMBOLS

| | |
|--|-----------------------------------|
| C. albicans Disks | C. albicans Disks |
| C. albicans 0.3% Sodium Hydroxide | C. albicans 0.3% Sodium Hydroxide |
| C. albicans Color Developer | C. albicans Color Developer |
| Plastic Tubes with Caps | Plastic Tubes with Caps |
| Applicator Sticks | Applicator Sticks |

6. PRECAUTIONS

- This product is for *in vitro* diagnostic use and should be used by properly trained individuals.
- Standard precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and test materials after use.

- Directions should be read and followed carefully.
- Caution! The Color Developer is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
- Refer to the Safety Data Sheet for detailed information on reagent chemicals.
- Reagents are provided at the necessary working strength and should be dispensed directly from the dropper bottles. Do not dilute reagents.
- Do not interchange reagents between kits of different lots.

Contents

Dimethyl sulfoxide 67-68-5
Acetic acid 64-19-7
2-Methoxyethanol 109-86-4
Sodium hydroxide 1310-73-2

DANGER

US ONLY

US & EU

| | |
|----------------|---|
| H315 | Causes skin irritation |
| H319 | Causes serious eye irritation |
| H335 | May cause respiratory irritation |
| H336 | May cause drowsiness or dizziness |
| H360 | May damage fertility. May damage the unborn child |
| H373 | May cause damage to organs through prolonged or repeated exposure |
| P201 | Obtain special instructions before use |
| P202 | Do not handle until all safety precautions have been read and understood |
| P281 | Use personal protective equipment as required |
| P264 | Wash face, hands and any exposed skin thoroughly after handling |
| P280 | Wear eye/face protection |
| P260 | Do not breathe dust/fume/gas/mist/vapors/spray |
| P271 | Use only outdoors or in a well-ventilated area |
| P308+P313 | IF exposed or concerned: Get medical attention/advice |
| P304+P340 | IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing |
| P302+P352 | IF ON SKIN: Wash with plenty of soap and water |
| P332+P313 | If skin irritation occurs: Get medical advice/attention |
| P362 | Take off contaminated clothing and wash before reuse |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing |
| P337+P313 | If eye irritation persists: Get medical advice/attention |
| P405 | Store locked up |
| P403+P233 | Store in a well-ventilated place. Keep container tightly closed |
| P501 | Dispose of contents/container to an approved waste disposal plant |

Hazards not otherwise classified (HNOC)- None identified
WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

7. STORAGE

Store product in the dark in its original container at 2-8°C until used. Do not freeze or overheat. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

8. PRODUCT DETERIORATION

This product should not be used if (1) the color of the disks has changed from white, (2) the expiration date has passed, or (3) there are other signs of deterioration. Protect disks from moisture by removing from the vial

only those disks necessary for testing. Promptly replace the cap and return the vial to 2-8°C.

Note: Precipitate formation or color variation from straw to pink in the Color Developer is common and does not affect test performance.

9. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.⁵

10. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Forceps, (7) Demineralized water, (8) Micropipette (5-10 μ l) or microbiological loop (10 μ l).

11. PROCEDURE

Test isolates should be less than 96 hours old and growing on commonly used media, such as blood agar, Sabouraud dextrose agar, or mycobiotic agar. Ensure the test isolate is in pure culture or select colonies that are well-isolated from surrounding organisms.

- Using forceps remove a Candida Albicans Test Kit disk from the vial and place it in one of the tubes provided.
- Add one drop of demineralized water to moisten the disk (~10 μ l). Do not oversaturate.
- With an applicator stick (provided) remove a heavy visible paste of the test isolate and rub it into the disk.
- Cap the tube and incubate at 35-37°C for 30 minutes.
- After incubation, add one drop of Sodium Hydroxide, 0.3% and examine tube against a white background for a distinct yellow color. Record the presence or absence of color.
- Add one drop of the Color Developer and examine for a pink to red color within one minute. Record the presence or absence of color.

12. INTERPRETATION

NGL:

Positive Test - Distinct yellow color after addition of Sodium Hydroxide, 0.3%

Negative Test - No yellow color after addition of Sodium Hydroxide, 0.3%

PRO:

Positive Test - Pink to red color within one minute of adding Color Developer

Negative Test - No color change or slight yellow color within one minute of adding Color Developer

13. EXPECTED VALUES

C. albicans is positive for both tests. If either the NGL or the PRO test is negative, or if both tests are negative, further testing is required for identification of the yeast isolate.

Note: *Candida dubliniensis* is also positive for both NGL and PRO. Growth at elevated temperatures and colony morphology on differential media have been reported to facilitate differentiation of *C. albicans* and *C. dubliniensis*.^{6,7}

14. QUALITY CONTROL

All lot numbers of Candida Albicans Test Kit have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

| CONTROL | INCUBATION | RESULTS | |
|--|-------------------------------|---------|-----|
| | | NGL | PRO |
| <i>Candida Albicans</i> ATCC® 60193 | Ambient, 30 min. @ 35-37°C | + | + |
| <i>Candida krusei</i> ATCC® 14243 | Ambient, 30 min. @ 35-37°C | — | — |

IFU30851401, Revised March, 2021 Printed in UK

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS
66215, USA

www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939

www.oxid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

15. LIMITATIONS

1. False-negative reactions may occur with *Candida Albicans* Test Kit if inadequate inoculum is used.
2. *C. dubliniensis*, first described in 1995, also produces germ tubes and chlamydospores.⁸ Growth at elevated temperatures and morphology on differential media have been shown to facilitate differentiation of *C. albicans* and *C. dubliniensis*.^{6,7}
3. This test is only part of the overall scheme for identification of yeasts. Further testing is required for definitive identification. Consult appropriate references for further instructions.^{1,5}

16. PERFORMANCE CHARACTERISTICS

In a comparative study, *Candida Albicans* Test Kit was evaluated in parallel with conventional methods and a commercially available identification system. A total of 583 fresh clinical isolates were tested, including *C. albicans* (303 isolates) and 10 other species of yeasts. The study results demonstrated *Candida Albicans* Test Kit has a sensitivity of 98.7% and a specificity of 99.6%.⁹

17. BIBLIOGRAPHY

1. Winn, W.G., S.D. Allen, J. William, E.W. Koneman, G. Procop, P.C. Schreckenberger, and G. Woods. 2006. *Koneman's Color Atlas and Textbook of Diagnostic Microbiology*. 6th ed. Lippincott Williams and Wilkins, Baltimore, MD.
2. Perry, J.L. and G.R. Miller. 1987. *J. Clin. Microbiol.* 25:2424-2425.
3. Perry, J.L., G.R. Miller, and D.L. Carr. 1990. *J. Clin. Microbiol.* 28:614-615.
4. Dealler, S.F. 1991. *J. Clin. Microbiol.* 29:1081-1082.
5. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller. 2007. *Manual of Clinical Microbiology*. 9th ed. ASM Press, Washington, D.C.
6. Al Mosaid, A., D. Sullivan, I.F. Salkin, D. Shanley, and D.C. Coleman. 2001. *J. Clin. Microbiol.* 39:323-327.
7. Wabale, V.R., A.S. Kagal, R.S. Mani, and R. Bharadwaj. 2007. *Indian J. Med. Microbiol.* 25:304-305. Retrieved October 1, 2008 from: <http://www.ijmm.org/text.asp?2007/25/3/304/34787>.
8. Sullivan, D.J., T.J. Westerneng, K.A. Haynes, D.E. Bennett, and D.C. Coleman. 1995. 141:15071521.
9. Crist, A.E., Jr., T.J. Dietz, and K. Kampschroer. 1996. *J. Clin. Microbiol.* 34:2616-2618.

18. PACKAGING

REF R30851401, *Candida Albicans* Test Kit ..50 Tests/Kit

19. SYMBOL LEGEND

| | |
|------------|---|
| REF | Catalogue Number |
| IVD | <i>In Vitro</i> Diagnostic Medical Device |
| | Consult Instructions for Use (IFU) |
| | Temperature Limitations (Storage temp.) |
| LAB | For Laboratory Use Only |
| LOT | Batch Code (Lot Number) |
| | Use By (Expiration Date) |
| | Manufactured by |

ATCC is a registered trademark of American Type Culture Collection.

For technical information contact your local distributor.

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA