# **Dehydrated Culture Media**

MRS BROTH (DE MAN,

You are viewing the printer friendly version of this page. To return to the regular view click here.

### ROGOSA, SHARPE)

Code: CM0359

A non-selective medium for profuse growth of `lactic acid bacteria'.

#### Typical Formula\*

gm/litre

Peptone

10.0

'Lab-Lemco' powder

Yeast extract

1

Glucose

20.0

Sorbitan mono-oleate

1 ml

Dipotassium hydrogen phosphate

2.0

Sodium acetate 3H<sub>2</sub>O

5.0

Triammonium citrate

2.0

Magnesium sulphate 7H2O

0.2

Manganese sulphate 4H<sub>2</sub>O

0.05

pH 6.2 ± 0.2 @ 25°C

#### Directions

Add 52g to 1 litre of distilled water at approximately 60°C. Mix until completely dissolved. Dispense into final containers and sterilise by autoclaving at 121°C for 15 minutes.

#### Description

MRS Broth may be used for tests in the identification of lactobacilli, such as temperature dependence, growth in 4% NaCl, growth in 0.4% Teepol, etc. as recommended by Sharpe, Fryer and Smith<sup>1</sup>.

#### Storage conditions and Shelf life

Store the dehydrated medium at 10-30°C and use before the expiry date on the label. Store the prepared medium at 2-8°C.

#### **Appearance**

Dehydrated medium: Dark straw powder Prepared medium: Amber coloured solution

## **Quality control**

Positive control: Expected results
Lactobacillus gasseri ATCC® 19992 \* Turbid growth

Negative control:

Uninoculated medium No change

#### Reference

1. Sharpe M. Elisabeth, Fryer T. F. and Smith D. G. (1966) `Identification of the Lactic Acid Bacteria' in `Identification Method for Microbiologists Part A' (Gibbs B. M. and Skinner F. A. eds.) London and New York, Academic Press. Pages 65-79.

©2001 - 2020 Oxoid Limited, All rights reserved.
Copyright, Disclaimer and Privacy Policy | Conditions of Sale | About Us | Cookies Thermo Fisher Scientific Inc.

<sup>\*</sup> Adjusted as required to meet performance standards

<sup>\*</sup> This organism is available as a Culti-Loop®