

THE EXPERTISE OF UPSTATE®, CHEMICON® AND LINCO®
IS NOW A PART OF MILLIPORE



## RABBIT ANTI-HUMAN INTEGRIN alpha2 [CD49b] POLYCLONAL ANTIBODY

**CATALOG NUMBER:** AB1936 QUANTITY: 100 μL

LOT NUMBER: 2139993

**ALTERNATE NAMES:** CD49b **EPITOPE:** Cytoplasmic domain

Human alpha2 integrin subunit. Specificity determined by immunoprecipitation from a SPECIFICITY:

detergent extract of <sup>35</sup>S methionine labeled endothelial cells. The epitope recognized is on the cytoplasmic domain of the protein and is thus intracellular. Shows no cross-reactivity to

alpha1, alpha3, alpha4, alpha5, alpha6, or alphaV subunits.

**IMMUNOGEN:** Synthetic peptide derived from the C-terminal sequence of the human alpha2 integrin

subunit.

**APPLICATIONS:** Immunofluorescence: 1:500 (suggested for use on acetone-fixed tissue only)

ELISA/RIA: 1:1000

Western blot: 1:1000 The antibody has been used to report α2 expression in HUVECs

during tubule-like structure formation in collagen (Deroanne, 1996). Immunoprecipitation: 5  $\mu$ L precipitates  $\alpha 2/\beta 1$  from 5 x 10<sup>6</sup> cells.

Immunohistochemistry: 1:500

Optimal working dilutions must be determined by end user.

SPECIES REACTIVITY: Human, rat, mouse, hamster, and chicken. Reactivity with other species has not been

determined.

Neat rabbit antiserum FORMAT:

PRESENTATION: Liquid with 0.05% sodium azide

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 12 months upon date of receipt. Avoid

repeated freeze/thaw cycles.

**REFERENCES:** Plantefaber, LC, et al., Cell (1989) 56:281-290.

> J Biol Chem, (1992) 267:18303-18307. Mol Biol of the Cell, (1993) 4:973-982.

Deroanne, CF, et al., Exper Cell Res, (1996) 224:215-223. Valles, AM, et al., Cell Adhes Commun, (1996) 4:187-199. Wang, XQ, et al., Mol Biol of the Cell, (1998) 9: 865-874.

Sherry, DM, et al, Graefe's Arch Clin Exp Ophthalmol, (2001) 239:278-287.





Important Note:

During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2012 Merck KGaA, Darmstadt. All rights reserved. No part of these works may be reproduced in any form without permission in writing.