

Anti-Sox9

Polyclonal Antibody



Cat. # AB5535

Lot # 3063352

FOR RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC PROCEDURES
NOT FOR HUMAN OR ANIMAL CONSUMPTION

Pack Size: 100 µg

Concentration 1.0 mg/ml

Storage: 2-8°C

Certificate of Analysis

Your Partner in Research
1-800-MILLIPORE
www.emdmillipore.com

page 1 of 3

Applications	Species Cross-Reactivity	Antibody Isotype	Epitope/Region	Host Species	Molecular Weight	Accession #
WB, IH, ChIP, ChIP-seq, ICC, IF	H, M, R, Ch, (Bo), (Eq), (Fe), (Ov)	N/A	C-terminal	Rb	~65 kDa	P48436

Immunogen	KLH-conjugated linear peptide corresponding to the C-terminal sequence of human Sox9.
Specificity	Recognizes Sox9.
Species Cross-reactivity	Human (H), Mouse (M), Rat (R), Chicken (Ch). Predicted to react with Bovine (Bo), Equine (Eq), Feline (Fe), and Ovine (Ov) based on 100% sequence homology.
Molecular Weight	~56-65 kDa observed.
Method of Purification	Affinity purified.
Presentation	Purified rabbit polyclonal antibody in buffer containing 0.1M Tris-Glycine (pH 7.4) 150mM NaCl with 0.05% sodium azide without glycerol.
Storage and Handling	Stable for 6 months at 2-8°C in undiluted aliquots from date of receipt.
Quality Control Testing	Evaluated by Western Blotting in L6 cell lysate. Western Blotting Analysis (WB): An 1:2000 dilution of this antibody detected Sox9 in L6 rat skeletal muscle myoblast lysate.
Additional Applications	Immunohistochemistry (IH) Analysis: An 1:1,000 dilution from a representative lot detected Sox9 in murine embryonic bone and adult cartilage tissue sections.

(Continued on page 2)

upstate | CHEMICON | Linco

© 2011 EMD Millipore Corporation, a division of Merck KGaA, Darmstadt, Germany. All rights reserved. Millipore, the M mark, Upstate, Chemicon, Linco and all other trademarks, unless specifically identified above in the text as belonging to a third party, are registered trademarks of Merck KGaA, Darmstadt, Germany.

Submit your published journal article, and earn credit toward future EMD Millipore purchases. Visit www.emdmillipore.com/publicationrewards to learn more!

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited. Millipore, the M mark, Upstate, Chemicon, Linco and all other trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2018 Merck KGaA, Darmstadt, Germany. All rights reserved.

Rev.6.0 /2017-03-02/AB5535CA/MA

Additional Applications

Western Blotting Analysis: An 1:500 dilution from a representative lot detected Sox9 in human PC3 prostate cancer cells and HepG2 hepatocytes.

Chromatin Immunoprecipitation (ChIP) Analysis: A representative lot detected Sox9 occupancy at target chromatin sites by ChIP using chromatin preparations from P1 post-natal mouse rib chondrocytes (Ohba, S., et al. (2015). *Cell Rep.* 12(2):229-243).

Chromatin Immunoprecipitation (ChIP) Analysis: A representative lot detected Sox9 occupancy at the Bmi promoter in Z/sox9tg but not in wild-type control mouse embryonic fibroblasts/MEFs (Matheu, A., et al. (2012). *Cancer Res.* 72(5):1301-1315).

ChIP-sequencing (ChIP-seq) Analysis: A representative lot detected Sox9-targeted chromatin sites by a genome-wide ChIP-seq analysis using chromatin preparations from P1 post-natal mouse rib chondrocytes (Ohba, S., et al. (2015). *Cell Rep.* 12(2):229-243)

Immunofluorescence (IF) Analysis: A representative lot detected the accumulation of Sox9-positive oval cells by fluorescent immunohistochemistry staining of paraffin-embedded liver sections from transgenic mice treated with diethylnitrosamine to induce conditional liver HNF4a knockout (Saha, S.K., et al. (2014). *Nature.* 513(7516):110-114).

Immunofluorescence Analysis: Representative lots detected Sox9 immunoreactivity in paraffin-embedded mouse embryo sections by fluorescent immunohistochemistry (Carrasco, M., et al. (2012). *J. Clin. Invest.* 122(10):3504-3515; Sylva, M., et al. (2011). *PLoS One.* 6(8):e22616).

Immunofluorescence Analysis: Representative lots immunostained Müller glial cells in frozen mouse and chicken retina sections by fluorescent immunohistochemistry staining of (Muranishi, Y., and Furukawa, T. (2012). *J. Biomed. Biotechnol.* 2012:973140; Fischer, A.J., et al. (2011). *Neuroscience.* 178:250-260).

Immunohistochemistry Analysis: A representative lot immunostained the supporting cells (Sertoli) of the seminiferous tubules by immunohistochemistry staining of paraffin-embedded mouse testis sections (O'Shaughnessy, P.J., et al. (2012). *PLoS One.* 7(4):e35136).

(Continued on page 3)

Please visit www.emdmillipore.com for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

Millipore, the M mark, Upstate, Chemicon, Linco and all other trademarks, unless specifically identified above in the text as

belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2018 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified
Renewable Energy

Additional Applications

Immunocytochemistry Analysis: A representative lot detected the stem cell marker Sox9 by fluorescent immunocytochemistry staining of paraformaldehyde-fixed E-Cad/Lgr6+ human lung stem cells (HLSCs) clonally derived and passaged in culture (Oeztuerk-Winder, F., et al. (2012). EMBO J. 31(16):3431-3441).

Immunohistochemistry Analysis: A representative lot detected Sox9 immunoreactivity in various formalin-fixed, paraffin-embedded human tumor tissue sections (Matheu, A., et al. (2012). Cancer Res. 72(5):1301-1315).

Western Blotting Analysis: A representative lot detected the stem cell marker Sox9 in E-Cad/Lgr6+ human lung stem cells (HLSCs) clonally derived and passaged in culture (Oeztuerk-Winder, F., et al. (2012). EMBO J. 31(16):3431-3441).

Western Blotting Analysis: A representative lot detected upregulated Sox9 expression level in human colorectal cancer cell lines, HCT116, DLD1, and SW620 (Matheu, A., et al. (2012). Cancer Res. 72(5):1301-1315).

For sample data please visit – www.emdmillipore.com

Please visit www.emdmillipore.com for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

Millipore, the M mark, Upstate, Chemicon, Linco and all other trademarks, unless specifically identified above in the text as

belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2018 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified
Renewable Energy