

PE anti-human CD49d Antibody

304303 / 25 tests Catalog# / Size

304304 / 100 tests

9F10 Clone Workshop V S215

Other Names VLA-4 α chain, α4 integrin, Integrin α4 chain, ITGA4

Mouse IgG1, κ Isotype

Description CD49d is a 150 kD α integrin chain known as α_4 integrin or VLA-4 α chain. It forms a

heterodimer with either integrin $\beta1$ ($\alpha_4\beta_1$, VLA-4) or $\beta7$ ($\alpha_4\beta_7$). CD49d is expressed broadly on T lymphocytes, B lymphocytes, monocytes, thymocytes, eosinophils, basophils, mast cells, NK cells, dendritic cells, and some non-hematopoietic cells, but not on normal red blood cells, platelets or neutrophils. VLA-4 binds to VCAM-1 (CD106) and fibronectin. $\alpha_4\beta_7$ is the receptor for VCAM-1 and MAdCAM-1. CD49d participates in mononuclear cell trafficking to endothelial sites of inflammation and has roles in cell-cell interactions and cell adhesion to extracellular matrices. CD49d is involved in lymphocyte migration, T cell activation, and

hematopoietic stem cell differentiation. CD49d is a marker to isolate pure populations of Treg

cells due to its absence on Foxp3+ cells.

Product Details

Human, African Green, Baboon, Cat (Feline), Cattle (Bovine, Cow), Chimpanzee, Common Marmoset, Reactivity

Cynomolgus, Dog (Canine), Horse (Equine), Rhesus, Sheep (Ovine), Squirrel Monkey

Monoclonal Antibody Type

Host Species Mouse

Formulation Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA (origin

Preparation The antibody was purified by affinity chromatography, and conjugated with PE under optimal

conditions.

Concentration Lot-specific (please contact technical support for concentration and total µg amount, or use our Lookup

tool if you have a lot number.)

Storage & Handling The antibody solution should be stored undiluted between 2°C and 8°C, and protected from

prolonged exposure to light. Do not freeze.

Application FC - Quality tested

Recommended Usage Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric

analysis. For flow cytometric staining, the suggested use of this reagent is 5 μ l per million cells in 100

μl staining volume or 5 μl per 100 μl of whole blood.

Excitation Laser Blue Laser (488 nm)

Green Laser (532 nm)/Yellow-Green Laser (561 nm)

Application Notes Additional reported applications (for the relevant formats) include: immunohistochemical staining of

acetone-fixed frozen tissue sections, and in vitro T cell costimulation^{2,3}. The LEAF™ Purified antibody (Endotoxin <0.1 EU/ μ g, Azide-Free, 0.2 μ m filtered) is recommended for functional assays (Cat. No.

304310).

Application References

1. Schlossman S, et al. Eds. 1995. Leucocyte Typing V. Oxford University Press. New York.

(PubMed link indicates **BioLegend citation)**

- Jeong SH, et al. 2004. J. Virol. 78:6995. (Costim) Vogel TU, et al. 2002. J. Immunol. 169:4511. (Costim)
- Kleinewietfeld M, et al. 2009. Blood 113:827. (FC) PubMed
- Palacious F, et al. 2010. Blood 115:4488. PubMed
- 6. Yoshino N, et al. 2000. Exp. Anim. (Tokyo) 49:97. (FC)
- Sestak K, et al. 2007. Vet. Immunol. Immunopathol. 119:21.

Mattapallil MJ, et al. 2011. J. Immunol. 187:1977. PubMed

Product Citations

- 1. Ito K. et al. 2008. Matrix Biol. 28:11. PubMed
- Zhang Z, Shively J 2013. PLoS One. 8:67649. PubMed
- 3. Velázquez-Avila M. et al. 2019, Leukemia, 33:1337, PubMed

Antigen Details

Structure Integrin, type I transmembrane glycoprotein, 150 kD.

Distribution T cells, B cells, NK , dendritic cells, thymocytes, monocytes, eosinophils, mast cells.

Function Lymphocyte migration, T cell activation, stem cell differentiation.

Ligand/Receptor Fibronectin, VCAM-1, MAdCAM-1

Cell Type B cells, Dendritic cells, Eosinophils, Mast cells, Monocytes, NK cells, T cells, Thymocytes, Tregs

Biology Area Cell Adhesion, Cell Biology, Immunology, Innate Immunity

Molecular Family Adhesion Molecules, CD Molecules

Antigen References 1. Elices M, Ed.1995. Springer Semin. Immunopathol. 16(4).

2. Lobb RR and Helmer ME. et al. 1994. J. Clin. Invest. 94:1722.

Gene ID 3676


Related Protocols

Cell Surface Flow Cytometry Staining Protocol

Other Formats

APC anti-human CD49d, PE/Cyanine5 anti-human CD49d, Purified anti-human CD49d, Alexa Fluor® 594 anti-human CD49d, PerCP/Cyanine5.5 anti-human CD49d, PE/Cyanine7 anti-human CD49d, FITC anti-human CD49d, Brilliant Violet 510™ anti-human CD49d, Brilliant Violet 421™ anti-human CD49d, Per/Dazzle™ 594 anti-human CD49d, APC/Cyanine7 anti-human CD49d, Brilliant Violet 711™ anti-human CD49d, APC/Cyanine7 anti-human CD49d, Brilliant Violet 711™ anti-human CD49d, Alexa Fluor® 647 anti-human CD49d, Biotin anti-human CD49d, TotalSeq™-A0576 anti-human CD49d, Brilliant Violet 785™ anti-human CD49d, APC/Fire™ 750 anti-human CD49d, Ultra-LEAF™ Purified anti-human CD49d, TotalSeq™-C0576 anti-human CD49d

Product Data



Human peripheral blood lymphocytes stained with 9F10 PE

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587